



Acpr30 (49875): sc-73669

BACKGROUND

Acpr30 (adipocyte complement-related protein of 30 kDa or AdipoQ) is a secretory protein made exclusively in adipocytes with mRNA induced over 100-fold during adipocyte differentiation. Posttranscriptional modification of Acpr30 yields several oligomeric forms of varying molecular weight, including a monomer, a dimer, a trimer, a hexamer and a polymer. Acpr30 is an abundant serum protein, secreted exclusively from fat cells, and is implicated in energy homeostasis and obesity. Due to the dysregulation of Acpr30 in cases of obesity in humans and mice and the strong structural similarity to TNF α , Acpr30 is a suspected regulator of whole body energy homeostasis. In addition, regulated exocytosis of Acpr30 appears to require phosphatidylinositol-3-kinase activity, since insulin-stimulated Acpr30 secretion is blocked by pharmacologic inhibitors of this enzyme.

REFERENCES

1. Scherer, P.E., Williams, S., Fogliano, M., Baldini, G. and Lodish, H.F. 1995. A novel serum protein similar to C1q, produced exclusively in adipocytes. *J. Biol. Chem.* 270: 26746-26749.
2. Shapiro, L. and Scherer, P.E. 1998. The crystal structure of a complement-1q family protein suggests an evolutionary link to tumor necrosis factor. *Curr. Biol.* 8: 335-338.
3. Bogan, J.S. and Lodish, H.F. 1999. Two compartments for insulin-stimulated exocytosis in 3T3-L1 adipocytes defined by endogenous ACRP30 and Glut4. *J. Cell Biol.* 146: 609-620.
4. Kappes, A. and Loffler, G. 2000. Influences of ionomycin, dibutyryl-cyclo-AMP and tumour necrosis factor-alpha on intracellular amount and secretion of apM1 in differentiating primary human preadipocytes. *Horm. Metab. Res.* 32: 548-554.
5. Das, K., Lin, Y., Widen, E., Zhang, Y. and Scherer, P.E. 2001. Chromosomal localization, expression pattern, and promoter analysis of the mouse gene encoding adipocyte-specific secretory protein Acpr30. *Biochem. Biophys. Res. Commun.* 280: 1120-1129.
6. Fruebis, J., Tsao, T.S., Javorschi, S., Ebbets-Reed, D., Erickson, M.R., Yen, F.T., Bihain, B.E. and Lodish, H.F. 2001. Proteolytic cleavage product of 30 kDa adipocyte complement-related protein increases fatty acid oxidation in muscle and causes weight loss in mice. *Proc. Natl. Acad. Sci. USA* 98: 2005-2010.
7. Wang, Y., Lam, K.S., Chan, L., Chan, K.W., Lam, J.B., Lam, M.C., Hoo, R.C., Mak, W.W., Cooper, G.J. and Xu, A. 2006. Post-translational modifications of the four conserved lysine residues within the collagenous domain of adiponectin are required for the formation of its high molecular weight oligomeric complex. *J. Biol. Chem.* 281: 16391-16400.

CHROMOSOMAL LOCATION

Genetic locus: ADIPOQ (human) mapping to 3q27; Adipoq (mouse) mapping to 16.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

Acpr30 (49875) is a rat monoclonal antibody raised against full length recombinant Acpr30 of mouse origin.

PRODUCT

Each vial contains 100 μ g IgG2a in 1.0 ml PBS with < 0.1% sodium azide and protein stabilizer.

APPLICATIONS

Acpr30 (49875) is recommended for detection of Acpr30 of mouse and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Acpr30 siRNA (h): sc-43600, Acpr30 siRNA (m): sc-45891, Acpr30 shRNA Plasmid (h): sc-43600-SH, Acpr30 shRNA Plasmid (m): sc-45891-SH, Acpr30 shRNA (h) Lentiviral Particles: sc-43600-V and Acpr30 shRNA (m) Lentiviral Particles: sc-45891-V.

Molecular Weight of LMW Acpr30 monomer: 30 kDa.

Molecular Weight of HMW Acpr30 polymer: greater than 250 kDa.

Positive Controls: 3T3-L1 cell lysate: sc-2243 or MIA PaCa-2 cell lysate: sc-2285.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Immunofluorescence: use goat anti-rat IgG-FITC: sc-2011 (dilution range: 1:100-1:400) or goat anti-rat IgG-TR: sc-2782 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.