

CapZ- α (P-16): sc-74360

BACKGROUND

The F-actin family of capping proteins includes CapZ- α 1, CapZ- α 2, CapZ- α 3 and CapZ- β 3, all of which function in a calcium-dependent manner and bind to the fast growing barbed end of actin filaments, thereby blocking protein exchange at these ends. The F-actin capping protein complex is a heterodimer consisting of α and β subunits that caps the barbed ends of actin filaments and nucleates the polymerization of Actin monomers, yet does not sever actin filaments. CapZ- α 1, also known as F-actin-capping protein subunit α 1, is a 286 amino acid subunit of the heterodimer that forms the F-actin capping protein complex. CapZ- α 1 also has been shown to bind S-100 β chain, a signaling molecule involved in the calcium-sensitive assembly of intermediate filaments that has been linked to Alzheimer's disease.

REFERENCES

- Casella, J.F. and Torres, M.A. 1994. Interaction of CapZ with actin. The NH₂-terminal domains of the α 1 and β subunits are not required for actin capping, and α 1 β and α 2 β heterodimers bind differentially to actin. *J. Biol. Chem.* 269: 6992-6998.
- Hart, M.C., et al. 1997. Vertebrates have conserved capping protein α isoforms with specific expression patterns. *Cell Motil. Cytoskeleton* 38: 120-132.
- Inman, K.G., et al. 2002. Solution NMR structure of S-100B bound to the high-affinity target peptide TRTK-12. *J. Mol. Biol.* 324: 1003-1014.
- Wear, M.A., et al. 2003. How capping protein binds the barbed end of the actin filament. *Curr. Biol.* 13: 1531-1537.
- Hutchings, N.J., et al. 2003. Linking the T cell surface protein CD2 to the actin-capping protein CAPZ via CMS and CIN85. *J. Biol. Chem.* 278: 22396-22403.
- Huang, S., et al. 2003. *Arabidopsis* capping protein (AtCP) is a heterodimer that regulates assembly at the barbed ends of actin filaments. *J. Biol. Chem.* 278: 44832-44842.
- Pyle, W.G. 2004. Searching for the missing link: a role for the actin capping protein in heart failure. *Can. J. Cardiol.* 20: 1429-1432.
- Katoh, M. 2005. Wnt-2B: comparative integromics and clinical applications (Review). *Int. J. Mol. Med.* 16: 1103-1108.
- Cooper, J.A. and Sept, D. 2008. New insights into mechanism and regulation of actin capping protein. *Int. Rev. Cell Mol. Biol.* 267: 183-206.

CHROMOSOMAL LOCATION

Genetic locus: CAPZA1 (human) mapping to 1p13.2, CAPZA2 (human) mapping to 7q31.2; Capza1 (mouse) mapping to 3 F2.2, Capza2 (mouse) mapping to 6 A2.

SOURCE

CapZ- α (P-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of CapZ- α 1 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74360 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CapZ- α (P-16) is recommended for detection of CapZ- α 1 and CapZ- α 2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CapZ- α (P-16) is also recommended for detection of CapZ- α 1 and CapZ- α 2 in additional species, including equine, canine, bovine, porcine and avian.

Molecular Weight of CapZ- α : 36 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **CapZ- α (H-9): sc-374302** or **CapZ- α (C-7): sc-376134**, our highly recommended monoclonal alternatives to CapZ- α (P-16).