

HADHSC (D-20): sc-74650

BACKGROUND

HADHSC (hydroxyacyl-coenzyme A (CoA) dehydrogenase, short chain), also known as HAD, HNF4, HADH1, SCHAD or M/SCHAD (medium and short chain L-3-hydroxyacyl-CoA dehydrogenase), is a mitochondrial matrix protein expressed in pancreas, liver, heart, kidney and skeletal muscle. HADHSC exists as a homodimer that participates in lipid metabolism and is essential for the β -oxidation of medium and short chain fatty acids. More specifically, HADHSC catalyzes the dehydrogenation of 3-hydroxyacyl-CoAs to their corresponding 3-ketoacyl-CoAs while NAD⁺ is simultaneously reduced to NADH. Defects in HADHSC can lead to HADH (3- α -hydroxyacyl-CoA dehydrogenase) deficiency and familial hyperinsulinemic hypoglycemia 4 (HNF4). HADH deficiency is characterized as a metabolic disorder with patients exhibiting hepatocerebralopathy, hypoglycemia, myopathy or cardiomyopathy and sometimes experiencing sudden death. HNF4 is a disorder characterized by elevated Insulin secretion that, if left untreated, can cause brain damage from recurrent hypoglycemia episodes.

CHROMOSOMAL LOCATION

Genetic locus: HADH (human) mapping to 4q25; Hadh (mouse) mapping to 3 G3.

SOURCE

HADHSC (D-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of HADHSC of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74650 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HADHSC (D-20) is recommended for detection of HADHSC of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HADHSC (D-20) is also recommended for detection of HADHSC in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for HADHSC siRNA (h): sc-75222, HADHSC siRNA (m): sc-75223, HADHSC shRNA Plasmid (h): sc-75222-SH, HADHSC shRNA Plasmid (m): sc-75223-SH, HADHSC shRNA (h) Lentiviral Particles: sc-75222-V and HADHSC shRNA (m) Lentiviral Particles: sc-75223-V.


Molecular Weight of HADHSC isoforms 1/2: 34/42 kDa.

Positive Controls: HADHSC (h2): 293T Lysate: sc-170845, Hela whole cell lysate: sc-2200 or Hep G2 cell lysate: sc-2227.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

DATA



HADHSC (D-20): sc-74650. Western blot analysis of HADHSC expression in non-transfected 293T: sc-117752 (A), human HADHSC transfected 293T: sc-170845 (B) and HeLa (C) whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **HADHSC (A-5): sc-376525** or **HADHSC (LB-7): sc-100472**, our highly recommended monoclonal alternatives to HADHSC (D-20).