

MPP3 (E-18): sc-74683

BACKGROUND

MPP3, also known as discs large homolog 3 (DLG3), is a 585 amino acid protein belonging to the MAGUK family. Proteins in the MAGUK family, including CASK, ZO-1, ZO-2, p55 and MPP2, are known to regulate cell proliferation, intracellular junctions and signaling pathways by interacting with the cytoskeleton. MPP3 contains one guanylate kinase-like domain, one PDZ domain, two L27 domains and one SH3 domain. Via the PDZ domain, MPP3 interacts with the C-terminus of TSLC1. MPP3 has also been found to interact with SR-2A and SR-4.

REFERENCES

1. Smith, S.A., Holik, P., Stevens, J., Mazoyer, S., Melis, R., Williams, B., White, R. and Albertsen, H. 1996. Isolation of a gene (DLG3) encoding a second member of the discs-large family on chromosome 17q12-q21. *Genomics* 31: 145-150.
2. Online Mendelian Inheritance in Man, OMIM™. 1999. Johns Hopkins University, Baltimore, MD. MIM Number: 601114: World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Kim, M., Jiang, L.H., Wilson, H.L., North, R.A. and Surprenant, A. 2001. Proteomic and functional evidence for a P2X7 receptor signalling complex. *EMBO J.* 20: 6347-6358.
4. Shin, B.K., Wang, H., Yim, A.M., Le Naour, F., Brichory, F., Jang, J.H., Zhao, R., Puravs, E., Tra, J., Michael, C.W., Misek, D.E. and Hanash, S.M. 2003. Global profiling of the cell surface proteome of cancer cells uncovers an abundance of proteins with chaperone function. *J. Biol. Chem.* 278: 7607-7616.
5. Fukuhara, H., Masuda, M., Yageta, M., Fukami, T., Kuramochi, M., Maruyama, T., Kitamura, T., Murakami, Y. and Masuda, M. 2003. Association of a lung tumor suppressor TSLC1 with MPP3, a human homologue of *Drosophila* tumor suppressor Dlg. *Oncogene* 22: 6160-6165.
6. Kantardzhieva, A., Alexeeva, S., Versteeg, I. and Wijnholds, J. 2006. MPP3 is recruited to the MPP5 protein scaffold at the retinal outer limiting membrane. *FEBS J.* 273: 1152-1165.
7. Carayol, N., Chen, J., Yang, F., Jin, T., Jin, L., States, D. and Wang, C.Y. 2006. A dominant function of IKK/NFκB signaling in global lipopolysaccharide-induced gene expression. *J. Biol. Chem.* 281: 31142-31151.
8. Gavarini, S., Becamel, C., Altier, C., Lory, P., Poncet, J., Wijnholds, J., Bockaert, J. and Marin, P. 2006. Opposite effects of PSD-95 and MPP3 PDZ proteins on serotonin 5-hydroxytryptamine_{2C} receptor desensitization and membrane stability. *Mol. Biol. Cell* 17: 4619-4631.
9. Heller, G., Fong, K.M., Girard, L., Seidl, S., End-Pfützenreuter, A., Lang, G., Gazdar, A.F., Minna, J.D., Zielinski, C.C. and Zöchbauer-Müller, S. 2006. Expression and methylation pattern of TSLC1 cascade genes in lung carcinomas. *Oncogene* 25: 959-968.

CHROMOSOMAL LOCATION

Genetic locus: MPP3 (human) mapping to 17q21.31; Mpp3 (mouse) mapping to 11 D.

SOURCE

MPP3 (E-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MPP3 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74683 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MPP3 (E-18) is recommended for detection of MPP3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MPP3 (E-18) is also recommended for detection of MPP3 in additional species, including equine, canine and porcine.

Suitable for use as control antibody for MPP3 siRNA (h): sc-75817, MPP3 siRNA (m): sc-75818, MPP3 shRNA Plasmid (h): sc-75817-SH, MPP3 shRNA Plasmid (m): sc-75818-SH, MPP3 shRNA (h) Lentiviral Particles: sc-75817-V and MPP3 shRNA (m) Lentiviral Particles: sc-75818-V.

Molecular Weight of MPP3: 70 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **MPP3 (D-8): sc-377522** or **MPP3 (PZ-2D21): sc-134390**, our highly recommended monoclonal alternatives to MPP3 (E-18).