

# NAP1L1 (T-14): sc-74738

## BACKGROUND

Proper nucleosome assembly is critical for compacting DNA into chromatin. NAP1 (nucleosome assembly protein 1) is a nuclear protein that acts as a transcriptional regulator and functions in nucleosome assembly. NAP1L1 (nucleosome assembly protein 1-like 1), also known as NRP, is a 391 amino acid member of the nucleosome assembly protein (NAP) family and may be involved in mediating chromatin formation. Localized to the nucleus and expressed throughout the body, NAP1L1 contains acidic domains which are thought to mediate NAP1L1-histone interaction. Due to its role in DNA replication, NAP1L1 is implicated as an important regulator of cell proliferation. NAP1L1 shares 54% sequence similarity with the *Saccharomyces cerevisiae* Nap1 protein and may be a genetic marker for intestinal carcinomas.

## REFERENCES

1. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 164060. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
2. Rehtanz, M., Schmidt, H.M., Warthorst, U. and Steger, G. 2004. Direct interaction between nucleosome assembly protein 1 and the papillomavirus E2 proteins involved in activation of transcription. *Mol. Cell Biol.* 24: 2153-2168.
3. Okuwaki, M., Kato, K., Shimahara, H., Tate, S. and Nagata, K. 2005. Assembly and disassembly of nucleosome core particles containing histone variants by human nucleosome assembly protein I. *Mol. Cell Biol.* 25: 10639-10651.
4. Kidd, M., Modlin, I.M., Mane, S.M., Camp, R.L., Eick, G. and Latich, I. 2006. The role of genetic markers—NAP1L1, MAGE-D2, and MTA1—in defining small-intestinal carcinoid neoplasia. *Ann. Surg. Oncol.* 13: 253-262.
5. Modlin, I.M., Kidd, M., Latich, I., Zikusoka, M.N., Eick, G.N., Mane, S.M. and Camp, R.L. 2006. Genetic differentiation of appendiceal tumor malignancy: a guide for the perplexed. *Ann. Surg.* 244: 52-60.
6. Al-Dhaheri, M.H., Shah, Y.M., Basrur, V., Pind, S. and Rowan, B.G. 2006. Identification of novel proteins induced by estradiol, 4-hydroxytamoxifen and acolbifene in T47D breast cancer cells. *Steroids* 71: 966-978.
7. Eckey, M., Hong, W., Papaioannou, M. and Baniahmad, A. 2007. The nucleosome assembly activity of NAP1 is enhanced by Alien. *Mol. Cell Biol.* 27: 3557-3568.

## CHROMOSOMAL LOCATION

Genetic locus: NAP1L1 (human) mapping to 12q21.2; Nap1l1 (mouse) mapping to 10 D1.

## SOURCE

NAP1L1 (T-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of NAP1L1 of human origin.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74738 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-74738 X, 200 µg/0.1 ml.

## APPLICATIONS

NAP1L1 (T-14) is recommended for detection of NAP1L1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NAP1L1 (T-14) is also recommended for detection of NAP1L1 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for NAP1L1 siRNA (h): sc-75871, NAP1L1 siRNA (m): sc-75872, NAP1L1 shRNA Plasmid (h): sc-75871-SH, NAP1L1 shRNA Plasmid (m): sc-75872-SH, NAP1L1 shRNA (h) Lentiviral Particles: sc-75871-V and NAP1L1 shRNA (m) Lentiviral Particles: sc-75872-V.

NAP1L1 (T-14) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of NAP1L1: 45 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.