

## PIST (G-19): sc-74836

### BACKGROUND

PIST (PDZ protein interacting specifically with TC10), also known as GOPC (golgi associated PDZ and coiled-coil motif containing), CAL or FIG, is a 462 amino acid protein that localizes to the cytoplasm, as well as to the membrane of the golgi apparatus and to the cell junction. Expressed ubiquitously and containing one PDZ (DHR) domain, PIST functions as a homooligomer that interacts with a variety of proteins and plays a role in intracellular protein trafficking and degradation. Additionally, PIST is thought to regulate ionic currents via membrane channel modification and may also play a role in autophagy. Chromosomal aberrations in the gene encoding PIST are found in glioblastoma multiform (GBM), a common and aggressive form of brain tumor, suggesting a role for mutated PIST in carcinogenesis. Three isoforms of PIST exist due to alternative splicing events.

### REFERENCES

1. Yao, R., Maeda, T., Takada, S. and Noda, T. 2001. Identification of a PDZ domain containing Golgi protein, GOPC, as an interaction partner of frizzled. *Biochem. Biophys. Res. Commun.* 286: 771-778.
2. Cheng, J., Wang, H. and Guggino, W.B. 2004. Modulation of mature cystic fibrosis transmembrane regulator protein by the PDZ domain protein CAL. *J. Biol. Chem.* 279: 1892-1898.
3. He, J., Bellini, M., Xu, J., Castleberry, A.M. and Hall, R.A. 2004. Interaction with cystic fibrosis transmembrane conductance regulator-associated ligand (CAL) inhibits beta1-adrenergic receptor surface expression. *J. Biol. Chem.* 279: 50190-50196.
4. Hicks, S.W. and Machamer, C.E. 2005. Isoform-specific interaction of Golgin-160 with the Golgi-associated protein PIST. *J. Biol. Chem.* 280: 28944-28951.
5. Wente, W., Stroh, T., Beaudet, A., Richter, D. and Kreienkamp, H.J. 2005. Interactions with PDZ domain proteins PIST/GOPC and PDZK1 regulate intracellular sorting of the somatostatin receptor subtype 5. *J. Biol. Chem.* 280: 32419-32425.
6. Ito, H., Iwamoto, I., Morishita, R., Nozawa, Y., Asano, T. and Nagata, K. 2006. Identification of a PDZ protein, PIST, as a binding partner for Rho effector Rhotekin: biochemical and cell-biological characterization of Rhotekin-PIST interaction. *Biochem. J.* 397: 389-398.
7. Li, X., Zhang, J., Cao, Z., Wu, J. and Shi, Y. 2006. Solution structure of GOPC PDZ domain and its interaction with the C-terminal motif of neuroigin. *Protein Sci.* 15: 2149-2158.
8. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 606845. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
9. Cushing, P.R., Fellows, A., Villone, D., Boisguerin, P. and Madden, D.R. 2008. The relative binding affinities of PDZ partners for CFTR: a biochemical basis for efficient endocytic recycling. *Biochemistry* 47: 10084-10098.

### CHROMOSOMAL LOCATION

Genetic locus: GOPC (human) mapping to 6q22.1; Gopc (mouse) mapping to 10 B3.

### SOURCE

PIST (G-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of PIST of human origin.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74836 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

PIST (G-19) is recommended for detection of PIST of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PIST (G-19) is also recommended for detection of PIST in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for PIST siRNA (h): sc-76148, PIST siRNA (m): sc-76149, PIST shRNA Plasmid (h): sc-76148-SH, PIST shRNA Plasmid (m): sc-76149-SH, PIST shRNA (h) Lentiviral Particles: sc-76148-V and PIST shRNA (m) Lentiviral Particles: sc-76149-V.

Molecular Weight of PIST isoforms: 60/46/35 kDa.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **PIST (A-7): sc-393026**, our highly recommended monoclonal alternative to PIST (G-19).