

GSTO1 (E-12): sc-74935

BACKGROUND

GSTO1 (glutathione S-transferase ω 1), also known as p28 or GSTTlp28, is a 241 amino acid protein that localizes to the cytoplasm and contains both an N-terminal and a C-terminal GST domain. Expressed ubiquitously with highest expression in heart, liver and skeletal muscle, GSTO1 exists as a homodimer that functions as both a Glutathione-dependent thiol transferase and a dehydroascorbate reductase. Specifically, GSTO1 catalyzes the reaction of Glutathione with a wide variety of organic compounds to form thioethers, a process that is essential for the metabolism and detoxification of a variety of xenobiotics and carcinogens. Human GSTO1 shares 70% sequence homology with its rodent counterpart, suggesting a conserved role between species. Polymorphisms in the gene encoding GSTO1 may be associated with the development of childhood acute lymphoblastic leukemia, Parkinson's disease and Alzheimer disease.

REFERENCES

1. Ishikawa, T., Casini, A.F. and Nishikimi, M. 1998. Molecular cloning and functional expression of rat liver Glutathione-dependent dehydroascorbate reductase. *J. Biol. Chem.* 273: 28708-28712.
2. Kodym, R., Calkins, P. and Story, M. 1999. The cloning and characterization of a new stress response protein. A mammalian member of a family of θ class glutathione S-transferase-like proteins. *J. Biol. Chem.* 274: 5131-5137.
3. Board, P.G., Coggan, M., Chelvanayagam, G., Easteal, S., Jermini, L.S., Schulte, G.K., Danley, D.E., Hoth, L.R., Griffor, M.C., Kamath, A.V., Rosner, M.H., Chrnyk, B.A., Perregaux, D.E., Gabel, C.A., Geoghegan, K.F. and Pandit, J. 2000. Identification, characterization, and crystal structure of the ω class Glutathione transferases. *J. Biol. Chem.* 275: 24798-24806.
4. Yin, Z.L., Dahlstrom, J.E., Le Couteur, D.G. and Board, P.G. 2001. Immunohistochemistry of ω class glutathione S-transferase in human tissues. *J. Histochem. Cytochem.* 49: 983-987.
5. Li, Y.J., Oliveira, S.A., Xu, P., Martin, E.R., Stenger, J.E., Scherzer, C.R., Hauser, M.A., Scott, W.K., Small, G.W., Nance, M.A., Watts, R.L., Hubble, J.P., Koller, W.C., Pahwa, R., Stern, M.B., Hiner, B.C., Jankovic, J., Goetz, C.G., Mastaglia, F., Middleton, L.T., Roses, A.D., Saunders, A.M., et al. 2003. Glutathione S-transferase ω -1 modifies age-at-onset of Alzheimer disease and Parkinson disease. *Hum. Mol. Genet.* 12: 3259-3267.
6. Whitbread, A.K., Tetlow, N., Eyre, H.J., Sutherland, G.R. and Board, P.G. 2003. Characterization of the human ω class Glutathione transferase genes and associated polymorphisms. *Pharmacogenetics* 13: 131-144.
7. Whitbread, A.K., Mellick, G.D., Silburn, P.A., Le Couteur, D.G. and Board, P.G. 2004. Glutathione transferase ω class polymorphisms in Parkinson disease. *Neurology* 62: 1910-1911.
8. Wahner, A.D., Glatt, C.E., Bronstein, J.M. and Ritz, B. 2007. Glutathione S-transferase μ , ω , π , and θ class variants and smoking in Parkinson's disease. *Neurosci. Lett.* 413: 274-278.
9. Pongstaporn, W., Pakakasama, S., Sanguansin, S., Hongeng, S. and Petmitr, S. 2008. Polymorphism of glutathione S-transferase ω gene: association with risk of childhood acute lymphoblastic leukemia. *J. Cancer Res. Clin. Oncol.* E-published.

CHROMOSOMAL LOCATION

Genetic locus: GSTO1 (human) mapping to 10q25.1.

SOURCE

GSTO1 (E-12) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of GSTO1 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-74935 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GSTO1 (E-12) is recommended for detection of GSTO1 of human and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for GSTO1 siRNA (h): sc-75207, GSTO1 shRNA Plasmid (h): sc-75207-SH and GSTO1 shRNA (h) Lentiviral Particles: sc-75207-V.

Molecular Weight of GSTO1: 31 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

MONOS
Satisfaction
Guaranteed

Try **GSTO1/2 (H-12): sc-166121**, our highly recommended monoclonal alternative to GSTO1 (E-12).