

Ran BP-17 siRNA (m): sc-76347

BACKGROUND

The transport of proteins and large RNAs through the nuclear pore complexes (NPC) is an energy-dependent and regulated process. The import of proteins with a nuclear localization signal (NLS) is accomplished by recognition of one or more clusters of basic amino acids by the importin- α/β complex. The small GTPase Ran plays a key role in NLS-dependent protein import. Ran BP-17 (Ran-binding protein 17) is a 1,088 amino acid protein that belongs to the importin- β superfamily of nuclear transport receptors. Ran BP-17 is highly expressed in testis, moderately expressed in pancreas and weakly expressed in heart, placenta, lung, liver, thyroid, spinal cord, trachea and adrenal gland. The Ran BP-17 protein binds to nucleoporins and the GTP-bound form of Ran. Human Ran BP-16 and Ran BP-17 share 66% amino acid sequence identity. The Ran BP-17 gene is conserved in chimpanzee, canine, bovine, mouse and rat, and maps to human chromosome 5q35.1.

REFERENCES

1. Saltman, D.L., et al. 1993. A physical map of 15 loci on human chromosome 5q23-q33 by two-color fluorescence *in situ* hybridization. *Genomics* 16: 726-732.
2. Koch, P., et al. 2000. Identification of a novel putative Ran-binding protein and its close homologue. *Biochem. Biophys. Res. Commun.* 278: 241-249.
3. Kutay, U., et al. 2000. Identification of two novel RanGTP-binding proteins belonging to the importin β superfamily. *J. Biol. Chem.* 275: 40163-40168.
4. Online Mendelian Inheritance in Man, OMIM[™]. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 606141. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Su, X.Y., et al. 2006. HOX11L2/TLX3 is transcriptionally activated through T-cell regulatory elements downstream of BCL11B as a result of the t(5;14)(q35;q32). *Blood* 108: 4198-4201.
6. Koolen, D.A., et al. 2006. Holoprosencephaly and preaxial polydactyly associated with a 1.24 Mb duplication encompassing FBXW11 at 5q35.1. *J. Hum. Genet.* 51: 721-726.

CHROMOSOMAL LOCATION

Genetic locus: Ranbp17 (mouse) mapping to 11 A4.

PRODUCT

Ran BP-17 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Ran BP-17 shRNA Plasmid (m): sc-76347-SH and Ran BP-17 shRNA (m) Lentiviral Particles: sc-76347-V as alternate gene silencing products.

For independent verification of Ran BP-17 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-76347A, sc-76347B and sc-76347C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Ran BP-17 siRNA (m) is recommended for the inhibition of Ran BP-17 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Ran BP-17 gene expression knockdown using RT-PCR Primer: Ran BP-17 (h)-PR: sc-76346-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.