



## VPS33B (S-20): sc-79353

### BACKGROUND

Vacuolar sorting proteins (VPSs) are required for proper trafficking of endocytic and biosynthetic proteins to the vacuole and play an important role in the budding process of cells. VPS33B (vacuolar protein sorting 33 homolog B) is a 617 amino acid protein that localizes to the cytoplasmic side of the peripheral membrane, as well as to the late endosomal membrane, and belongs to the STXBP/Sec1 family. Expressed ubiquitously with highest expression in testis and lowest expression in lung, VPS33B is thought to play a role in vesicle-mediated protein trafficking to lysosomal compartments and may also be involved in membrane docking events at late endosomes. Defects in the gene encoding VPS33B are the cause of arthrogryposis-renal dysfunction-cholestasis syndrome (ARC), an autosomal recessive disorder that is characterized by renal tubular dysfunction, neurogenic arthrogryposis multiplex congenita and neonatal cholestasis with bile duct hypoplasia.

### REFERENCES

1. Carim, L., Sumoy, L., Andreu, N., Estivill, X. and Escarceller, M. 2000. Cloning, mapping and expression analysis of VPS33B, the human orthologue of rat Vps33b. *Cytogenet. Cell Genet.* 89: 92-95.
2. Huizing, M., Didier, A., Walenta, J., Anikster, Y., Gahl, W.A. and Krämer, H. 2001. Molecular cloning and characterization of human VPS18, VPS11, VPS16, and VPS33. *Gene* 264: 241-247.
3. Gissen, P., Johnson, C.A., Morgan, N.V., Stapelbroek, J.M., Forshew, T., Cooper, W.N., McKiernan, P.J., Klomp, L.W., Morris, A.A., Wraith, J.E., McClean, P., Lynch, S.A., Thompson, R.J., Lo, B., Quarrell, O.W., Di Rocco, M., Trembath, R.C., Mandel, H., Wali, S., Karet, F.E., Knisely, A.S., Houwen, R.H., Kelly, D.A. and Maher, E.R. 2004. Mutations in VPS33B, encoding a regulator of SNARE-dependent membrane fusion, cause arthrogryposis-renal dysfunction-cholestasis (ARC) syndrome. *Nat. Genet.* 36: 400-404.
4. Lo, B., Li, L., Gissen, P., Christensen, H., McKiernan, P.J., Ye, C., Abdelhaleem, M., Hayes, J.A., Williams, M.D., Chitayat, D. and Kahr, W.H. 2005. Requirement of VPS33B, a member of the Sec1/Munc18 protein family, in megakaryocyte and platelet  $\alpha$ -granule biogenesis. *Blood* 106: 4159-4166.
5. Gissen, P., Johnson, C.A., Gentle, D., Hurst, L.D., Doherty, A.J., O'Kane, C.J., Kelly, D.A. and Maher, E.R. 2005. Comparative evolutionary analysis of VPS33 homologues: genetic and functional insights. *Hum. Mol. Genet.* 14: 1261-1270.
6. Hershkovitz, D., Mandel, H., Ishida-Yamamoto, A., Chefetz, I., Hino, B., Luder, A., Indelman, M., Bergman, R. and Sprecher, E. 2008. Defective lamellar granule secretion in arthrogryposis, renal dysfunction, and cholestasis syndrome caused by a mutation in VPS33B. *Arch Dermatol.* 144: 334-340.
7. Bach, H., Papavinasasundaram, K.G., Wong, D., Hmama, Z. and Av-Gay, Y. 2008. Mycobacterium tuberculosis virulence is mediated by PtpA dephosphorylation of human vacuolar protein sorting 33B. *Cell Host Microbe* 3: 316-322.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### CHROMOSOMAL LOCATION

Genetic locus: VPS33B (human) mapping to 15q26.1; Vps33b (mouse) mapping to 7 D3.

### SOURCE

VPS33B (S-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of VPS33B of human origin.

### PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-79353 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

### APPLICATIONS

VPS33B (S-20) is recommended for detection of VPS33B of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

VPS33B (S-20) is also recommended for detection of VPS33B in additional species, including equine, canine, bovine and avian.

Suitable for use as control antibody for VPS33B siRNA (h): sc-76905, VPS33B siRNA (m): sc-76906, VPS33B shRNA Plasmid (h): sc-76905-SH, VPS33B shRNA Plasmid (m): sc-76906-SH, VPS33B shRNA (h) Lentiviral Particles: sc-76905-V and VPS33B shRNA (m) Lentiviral Particles: sc-76906-V.

Molecular Weight of VPS33B: 71 kDa.

Positive Controls: Mouse testis extract: sc-2405 or rat testis extract: sc-2400.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.