

JAML (Y-15): sc-79550

BACKGROUND

JAML (adhesion molecule, interacts with CXADR antigen 1), also known as junctional adhesion molecule-like, dendritic cell-specific protein CREA7-1, UNQ722/PRO1387, CREA7-4, AMICA1, Gm638, CREA7-1, FLJ37080, MGC-118815 or MGC118814, is a 394 amino acid single-pass type I membrane protein and novel adhesion molecule with localization at the cell plasma membrane in regions of cell-to-cell contact but not free cell borders, implying that JAML participates in homophilic interactions. JAML belongs to the immunoglobulin superfamily and is expressed in peripheral blood leukocytes, bone marrow, spleen, lymph nodes and fetal thymus, liver and spleen. JAML assists in the migration of leukocytes through endothelial and epithelial tissues, mediates binding with CAR (coxsackie and adenovirus receptor) via its Ig-like V-type domain 2, and strengthens the binding of leukocytes to endothelial cells. Consisting of a transmembrane segment, cytoplasmic tail, 19-amino acid signal peptide and 2 extracellular immunoglobulin-like domains, JAML also undergoes alternative splicing resulting in 3 isoforms and is encoded by a gene mapping to human chromosome 11q23.3.

REFERENCES

1. Moog-Lutz, C., et al. 2003. JAML, a novel protein with characteristics of a junctional adhesion molecule, is induced during differentiation of myeloid leukemia cells. *Blood* 102: 3371-3378.
2. Zen, K., et al. 2005. Neutrophil migration across tight junctions is mediated by adhesive interactions between epithelial coxsackie and adenovirus receptor and a junctional adhesion molecule-like protein on neutrophils. *Mol. Biol. Cell* 16: 2694-2703.
3. Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 609770. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Mirza, M., et al. 2007. Coxsackie-virus and adenovirus receptor is upregulated in migratory germ cells during passage of the blood-testis barrier. *Endocrinology* 148: 5459-5469.
5. Luissint, A.C., et al. 2008. JAML-mediated leukocyte adhesion to endothelial cells is regulated in cis by $\alpha 4/\beta 1$ integrin activation. *J. Cell Biol.* 183: 1159-1173.

CHROMOSOMAL LOCATION

Genetic locus: AMICA1 (human) mapping to 11q23.3; Amica1 (mouse) mapping to 9 A5.2.

SOURCE

JAML (Y-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an extracellular domain of JAML of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-79550 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

JAML (Y-15) is recommended for detection of JAML of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

JAML (Y-15) is also recommended for detection of JAML in additional species, including canine.

Suitable for use as control antibody for JAML siRNA (h): sc-75353, JAML siRNA (m): sc-75354, JAML shRNA Plasmid (h): sc-75353-SH, JAML shRNA Plasmid (m): sc-75354-SH, JAML shRNA (h) Lentiviral Particles: sc-75353-V and JAML shRNA (m) Lentiviral Particles: sc-75354-V.


Molecular Weight of JAML: 60 kDa.

Positive Controls: I-11.15 whole cell lysate: sc-364370.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



JAML (Y-15): sc-79550. Western blot analysis of JAML expression in I-11.15 whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.