

TRIM14 (K-15): sc-79762

BACKGROUND

The tripartite motif (TRIM) family of proteins are characterized by a conserved TRIM domain that includes a coiled-coil region, a B-box type zinc finger, one RING finger and three zinc-binding domains. TRIM14 (tripartite motif-containing 14) is a 442 amino acid protein belonging to the TRIM/RBCC family. With highest expression in liver, TRIM14 contains one B box-type zinc finger and one B30.2/SPRY domain. Through its B-box zinc finger domain, TRIM14 inhibits the transcriptional activity of PU.1 in a dose-dependent manner. TRIM14 exists as two alternatively spliced isoforms that are designated isoform α and isoform β . TRIM14 is encoded by a gene located on human chromosome 9, which consists of about 145 million bases, 4% of the human genome and encodes nearly 900 genes. Chromosome 9 may play a role in gender determination and deletion of the distal portion of 9p can lead to development of male to female sex reversal, the phenotype of a female with a male X,Y genotype.

REFERENCES

- Borden, K.L. and Freemont, P.S. 1996. The RING finger domain: a recent example of a sequence-structure family. *Curr. Opin. Struct. Biol.* 6: 395-401.
- Reymond, A., Meroni, G., Fantozzi, A., Merla, G., Cairo, S., Luzi, L., Riganelli, D., Zanaria, E., Messali, S., Cainarca, S., Guffanti, A., Minucci, S., Pelicci, P.G. and Ballabio, A. 2001. The tripartite motif family identifies cell compartments. *EMBO J.* 20: 2140-2151.
- Shyu, H.W., Hsu, S.H., Hsieh-Li, H.M. and Li, H. 2001. A novel member of the RBCC family, Trif, expressed specifically in the spermatids of mouse testis. *Mech. Dev.* 108: 213-216.
- Shyu, H.W., Hsu, S.H., Hsieh-Li, H.M. and Li, H. 2003. Forced expression of RNF36 induces cell apoptosis. *Exp. Cell Res.* 287: 301-313.
- Meroni, G. and Diez-Roux, G. 2005. TRIM/RBCC, a novel class of "single protein RING finger" E3 ubiquitin ligases. *Bioessays* 27: 1147-1157.
- Sardiello, M., Cairo, S., Fontanella, B., Ballabio, A. and Meroni, G. 2008. Genomic analysis of the TRIM family reveals two groups of genes with distinct evolutionary properties. *BMC Evol. Biol.* 8: 225.
- Ozato, K., Shin, D.M., Chang, T.H. and Morse, H.C. 2008. TRIM family proteins and their emerging roles in innate immunity. *Nat. Rev. Immunol.* 8: 849-860.

CHROMOSOMAL LOCATION

Genetic locus: TRIM14 (human) mapping to 9q22.33; Trim14 (mouse) mapping to 4 B1.

SOURCE

TRIM14 (K-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of TRIM14 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-79762 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

TRIM14 (K-15) is recommended for detection of TRIM14 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for TRIM14 siRNA (h): sc-76736, TRIM14 siRNA (m): sc-76737, TRIM14 shRNA Plasmid (h): sc-76736-SH, TRIM14 shRNA Plasmid (m): sc-76737-SH, TRIM14 shRNA (h) Lentiviral Particles: sc-76736-V and TRIM14 shRNA (m) Lentiviral Particles: sc-76737-V.


Molecular Weight of TRIM14: 50 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DTA



TRIM14 (K-15): sc-79762. Western blot analysis of TRIM14 expression in HeLa whole cell lysate.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.