

BLAP75 (G-20): sc-79873

BACKGROUND

Bloom's syndrome (BS) is a rare human genetic disorder characterized by dwarfism, immunodeficiency, genomic instability and cancer predisposition. BS is a result of mutation in the BLM gene, which encodes a protein that forms a multienzyme complex with topoisomerase III α , replication protein A and BLAP75 (also designated RecQ-mediated genome instability protein 1 or RMI1). BLM maintains genome integrity and catalyzes Holliday-junction branch migration and the annealing of complementary single-stranded DNA molecules. BLAP75, an OB-fold nucleic acid binding domain, is essential for the stability of the BLM complex *in vivo*. Specifically, BLAP75 enhances the ability of the BLM-Topo III α pair to branch migrate the Holliday junction or dissolve the double Holliday junction structure to yield non-crossover recombinants. BLAP75 colocalizes with BLM in subnuclear foci in response to DNA damage, and its depletion impairs the recruitment of BLM to these foci.

REFERENCES

1. Cheok, C.F., Bachrati, C.Z., Chan, K.L., Ralf, C., Wu, L. and Hickson, I.D. 2005. Roles of the Bloom's syndrome helicase in the maintenance of genome stability. *Biochem. Soc. Trans.* 33: 1456-1459.
2. Yin, J., Sobock, A., Xu, C., Meetei, A.R., Hoatlin, M., Li, L. and Wang, W. 2005. BLAP75, an essential component of Bloom's syndrome protein complexes that maintain genome integrity. *EMBO J.* 24: 1465-1476.
3. Raynard, S., Bussen, W. and Sung, P. 2006. A double Holliday junction dissolvosome comprising BLM, topoisomerase III α , and BLAP75. *J. Biol. Chem.* 281: 13861-13864.
4. Wu, L., Bachrati, C.Z., Ou, J., Xu, C., Yin, J., Chang, M., Wang, W., Li, L., Brown, G.W. and Hickson, I.D. 2006. BLAP75/RMI1 promotes the BLM-dependent dissolution of homologous recombination intermediates. *Proc. Natl. Acad. Sci. USA* 103: 4068-4073.
5. Chan, K.L., North, P.S. and Hickson, I.D. 2007. BLM is required for faithful chromosome segregation and its localization defines a class of ultrafine anaphase bridges. *EMBO J.* 26: 3397-3409.
6. Bussen, W., Raynard, S., Busygina, V., Singh, A.K. and Sung, P. 2007. Holliday junction processing activity of the BLM-Topo III α -BLAP75 complex. *J. Biol. Chem.* 282: 31484-31492.
7. Singh, T.R., Ali, A.M., Busygina, V., Raynard, S., Fan, Q., Du, C.H., Andreassen, P.R., Sung, P. and Meetei, A.R. 2008. BLAP18/RMI2, a novel OB-fold-containing protein, is an essential component of the Bloom helicase-double Holliday junction dissolvosome. *Genes Dev.* 22: 2856-2868.
8. Raynard, S., Zhao, W., Bussen, W., Lu, L., Ding, Y.Y., Busygina, V., Meetei, A.R. and Sung, P. 2008. Functional role of BLAP75 in BLM-topoisomerase III α -dependent holliday junction processing. *J. Biol. Chem.* 283: 15701-15708.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: RMI1 (human) mapping to 9q21.32; Rmi1 (mouse) mapping to 13 B1.

SOURCE

BLAP75 (G-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of BLAP75 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-79873 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

Available as TransCruz reagent for Gel Supershift and ChIP applications, sc-79873 X, 200 μ g/0.1 ml.

APPLICATIONS

BLAP75 (G-20) is recommended for detection of BLAP75 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for BLAP75 siRNA (h): sc-72651, BLAP75 siRNA (m): sc-72652, BLAP75 shRNA Plasmid (h): sc-72651-SH, BLAP75 shRNA Plasmid (m): sc-72652-SH, BLAP75 shRNA (h) Lentiviral Particles: sc-72651-V and BLAP75 shRNA (m) Lentiviral Particles: sc-72652-V.

BLAP75 (G-20) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of BLAP75: 70 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.