

# CLK3 (Y-13): sc-79907

## BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine (Ser/Thr) protein kinases. CLK3 (Cdc-like kinase 3), also known as PHCLK3, is a 638 amino acid nuclear and cytoplasmic protein that belongs to the Ser/Thr protein kinase family. Functioning as a dual-specificity kinase, CLK3 catalyzes the ATP-dependent phosphorylation of arginine- and serine-rich (SR) splicing factor proteins, thereby regulating both their function and their intranuclear distribution. Via its enzymatic activity, CLK3 is thought to be one of several members of a network of regulatory proteins that control RNA splicing events. Four isoforms of CLK3 exist due to alternative splicing.

## REFERENCES

1. Becker, W., et al. 1996. cDNA cloning and characterization of rat Clk3, a LAMMER kinase predominately expressed in testis. *Biochim. Biophys. Acta* 1312: 63-67.
2. Duncan, P.I., et al. 1998. The CLK2 and CLK3 dual-specificity protein kinases regulate the intranuclear distribution of SR proteins and influence pre-mRNA splicing. *Exp. Cell Res.* 241: 300-308.
3. Menegay, H., et al. 1999. The dual specificity protein kinase CLK3 is abundantly expressed in mature mouse spermatozoa. *Exp. Cell Res.* 253: 463-473.
4. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602990. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. García-Sacristán, A., et al. 2005. Protein kinase CLK/STY is differentially regulated during erythroleukemia cell differentiation: a bias toward the skipped splice variant characterizes postcommitment stages. *Cell Res.* 15: 495-503.
6. Olsen, J.V., et al. 2006. Global, *in vivo*, and site-specific phosphorylation dynamics in signaling networks. *Cell* 127: 635-648.
7. Wissing, J., et al. 2007. Proteomics analysis of protein kinases by target class-selective prefractionation and tandem mass spectrometry. *Mol. Cell Proteomics* 6: 537-547.

## CHROMOSOMAL LOCATION

Genetic locus: CLK3 (human) mapping to 15q24.1; Clk3 (mouse) mapping to 9 B.

## SOURCE

CLK3 (Y-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of CLK3 of human origin.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-79907 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

CLK3 (Y-13) is recommended for detection of CLK3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CLK3 siRNA (h): sc-72925, CLK3 siRNA (m): sc-72926, CLK3 shRNA Plasmid (h): sc-72925-SH, CLK3 shRNA Plasmid (m): sc-72926-SH, CLK3 shRNA (h) Lentiviral Particles: sc-72925-V and CLK3 shRNA (m) Lentiviral Particles: sc-72926-V.

Molecular Weight of CLK3: 74 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker<sup>™</sup> compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz<sup>™</sup> Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.