



## DBP (2G2): sc-80624

### BACKGROUND

Vitamin D-binding protein (DBP) is a multi-functional serum protein that binds to the plasma membranes of numerous cell types and mediates a variety of cellular functions. The locus of the DBP protein (also known as group-specific component protein or GC) is located at human chromosome 4q13.3. DBP functions in organ-specific transportation of vitamin D and its metabolites to the various target organs of the vitamin D endocrine system. In addition, DBP has immunomodulatory properties and is able to bind to the surface of leukocytes. DBP binds to the plasma membrane through a chondroitin sulfate proteoglycan. DBP serves as a co-chemotactic factor for C5a to enhance the chemotactic activity of C5a. DBP can also bind to globular Actin with high affinity and is involved in the clearance of Actin from the blood. DBP plays an important role in osteoclast differentiation. The diverse cellular functions of DBP require its cell surface binding ability to mediate different biological processes.

### REFERENCES

- DiMartino, S.J., et al. 1999. Initial characterization of the vitamin D binding protein (Gc-globulin) binding site on the neutrophil plasma membrane: evidence for a chondroitin sulfate proteoglycan. *J. Immunol.* 163: 2135-2142.
- Pani, M.A., et al. 1999. Vitamin D binding protein alleles and susceptibility for type 1 diabetes in Germans. *Autoimmunity* 31: 67-72.
- Papiha, S.S., et al. 1999. Vitamin D binding protein gene in male osteoporosis: association of plasma DBP and bone mineral density with (TAAA)(n)-Alu polymorphism in DBP. *Calcif. Tissue Int.* 65: 262-266.
- Hirai, M., et al. 2000. Variations in vitamin D-binding protein (group-specific component protein) are associated with fasting plasma insulin levels in Japanese with normal glucose tolerance. *J. Clin. Endocrinol. Metab.* 85: 1951-1953.
- Swamy, N., et al. 2000. Probing with vitamin D sterol-binding pocket of human vitamin D-binding protein with bromoacetate affinity labeling reagents containing the affinity probe at C-1, C-6, C-11, and C-19 positions of parent vitamin D sterols. *Arch. Biochem. Biophys.* 373: 471-478.
- DiMartino, S.J., et al. 2001. Elastase controls the binding of the vitamin D-binding protein (Gc-globulin) to neutrophils: a potential role in the regulation of C5a co-chemotactic activity. *J. Immunol.* 166: 2688-2694.
- Bogaerts, I., et al. 2001. Purification, crystallization and preliminary X-ray investigation of the complex of human vitamin D binding protein and rabbit muscle Actin. *Acta Crystallogr. D Biol. Crystallogr.* 57: 740-742.
- Swamy, N., et al. 2001. Baculovirus-expressed vitamin D-binding protein-macrophage activating factor (DBP-maf) activates osteoclasts and binding of 25-hydroxyvitamin D(3) does not influence this activity. *J. Cell. Biochem.* 81: 535-546.

### CHROMOSOMAL LOCATION

Genetic locus: GC (human) mapping to 4q13.3.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### SOURCE

DBP (2G2) is a mouse monoclonal antibody raised against full length native DBP isolated from human plasma.

### PRODUCT

Each vial contains 100 µg IgG<sub>1</sub> in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

### APPLICATIONS

DBP (2G2) is recommended for detection of DBP of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DBP siRNA (h): sc-41375, DBP shRNA Plasmid (h): sc-41375-SH and DBP shRNA (h) Lentiviral Particles: sc-41375-V.

Molecular Weight of DBP: 58 kDa.

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotting A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.