

# Orexin-B (C-19): sc-8071

## BACKGROUND

The hypothalamus is essential for maintaining homeostasis by integrating the vertebrate endocrine and nervous systems, thereby controlling temperature, thirst and hunger. Orexin-A and Orexin-B (also designated hypocretins) are hypothalamic neuropeptides that are derived from a single precursor, prepro-orexin, by proteolytic processing. These peptides bind to and activate two closely related G protein-coupled receptors, designated Orexin receptor-1 and Orexin receptor-2. Orexin-A protein and prepro-orexin mRNA are localized to neurons within the lateral section of the hypothalamus, designated the "feeding center". Prepro-orexin mRNA is up-regulated during fasting, suggesting that orexins may play a role in the central feedback mechanism that regulates feeding behavior. Orexin has been shown to increase the release of GABA and Glutamate from axons, a response seen as a result of most synaptic activities in the hypothalamic region.

## CHROMOSOMAL LOCATION

Genetic locus: HCRT (human) mapping to 17q21.2; Hcrt (mouse) mapping to 11 D.

## SOURCE

Orexin-B (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of Orexin-B of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-8071 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Orexin-B (C-19) is recommended for detection of Orexin-B processed active peptide of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Orexin-B (C-19) is also recommended for detection of Orexin-B processed active peptide in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for Orexin-A/B siRNA (h): sc-42152, Orexin-A/B siRNA (m): sc-42153, Orexin-A/B shRNA Plasmid (h): sc-42152-SH, Orexin-A/B shRNA Plasmid (m): sc-42153-SH, Orexin-A/B shRNA (h) Lentiviral Particles: sc-42152-V and Orexin-A/B shRNA (m) Lentiviral Particles: sc-42153-V.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## SELECT PRODUCT CITATIONS

- Liu, R.J., et al. 2002. Hypocretins (Orexins) regulate serotonin neurons in the dorsal raphe nucleus by excitatory direct and inhibitory indirect actions. *J. Neurosci.* 22: 9453-9464.
- Downs, J.L., et al. 2007. Orexin neuronal changes in the locus coeruleus of the aging rhesus macaque. *Neurobiol. Aging* 28: 1286-1289.
- Suzuki, H., et al. 2007. Differential distribution of Orexin-A-like and Orexin R-1-like immunoreactivities in the *Xenopus* pituitary. *Tissue Cell* 39: 423-430.
- Su, J., et al. 2008. Distribution of Orexin-B and its relationship with GnRH in the pig hypothalamus. *Res. Vet. Sci.* 85: 315-823.
- Xie, X., et al. 2008. Hypocretin/Orexin and nociceptin/orphanin FQ coordinately regulate analgesia in a mouse model of stress-induced analgesia. *J. Clin. Invest.* 118: 2471-2481.
- Hsu, D.T., et al. 2009. Paraventricular thalamic nucleus: subcortical connections and innervation by serotonin, Orexin, and corticotropin-releasing hormone in macaque monkeys. *J. Comp. Neurol.* 512: 825-848.
- López, J.M., et al. 2009. Comparative immunohistochemical analysis of the distribution of orexins (hypocretins) in the brain of amphibians. *Peptides* 30: 873-887.
- López, J.M., et al. 2009. Distribution of orexin/hypocretin immunoreactivity in the brain of the lungfishes *Protopterus dolloi* and *Neoceratodus forsteri*. *Brain Behav. Evol.* 74: 302-322.
- Bullmann, T., et al. 2010. Expression of the embryonal isoform (ON/3R) of the microtubule-associated protein tau in the adult rat central nervous system. *J. Comp. Neurol.* 518: 2538-2553.
- Domínguez, L., et al. 2010. Immunohistochemical localization of orexins (hypocretins) in the brain of reptiles and its relation to monoaminergic systems. *J. Chem. Neuroanat.* 39: 20-34.
- Li, F.W., et al. 2011. Behavioural and neuronal activation after microinjections of AMPA and NMDA into the perifornical lateral hypothalamus in rats. *Behav. Brain Res.* 224: 376-386.
- Li, Y., et al. 2011. Orexins in the midline thalamus are involved in the expression of conditioned place aversion to morphine withdrawal. *Physiol. Behav.* 102: 42-50.
- Michinaga, S., et al. 2011. Orexin neurons in hypothalamic slice cultures are vulnerable to endoplasmic reticulum stress. *Neuroscience* 190: 289-300.


 MONOS  
 Satisfation  
 Guaranteed

Try **Orexin-A (KK09): sc-80263**, our highly recommended monoclonal alternative to Orexin-B (C-19).