

CYPOR (C-19): sc-8091

BACKGROUND

P450 enzymes constitute a family of monooxygenase enzymes that are involved in the metabolism of a wide array of endogenous and xenobiotic compounds. Several P450 enzymes have been classified by sequence similarities as members of the CYP1A and CYP2A subfamilies. CYPOR, also known as cytochrome P450 reductase and NADPH cytochrome P450 reductase, is a microsomal enzyme responsible for the transfer of electrons from NADPH to cytochrome P450 enzymes during the P450 catalytic cycle. CYPOR is localized to the endoplasmic reticulum, where it is also able to transfer electrons to heme oxygenase and cytochrome b5. CYPOR is structurally related to two separate flavoprotein families, ferredoxin nucleotide reductase (FNR) and flavodoxin. Electron transfer of CYPOR requires the binding of two flavin cofactors, FAD and FMN, to the FNR and flavodoxin domains, respectively.

CHROMOSOMAL LOCATION

Genetic locus: POR (human) mapping to 7q11.23; Por (mouse) mapping to 5 G2.

SOURCE

CYPOR (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of CYPOR of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-8091 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

CYPOR (C-19) is recommended for detection of CYPOR of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CYPOR (C-19) is also recommended for detection of CYPOR in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for CYPOR siRNA (h): sc-35147, CYPOR siRNA (m): sc-35148, CYPOR shRNA Plasmid (h): sc-35147-SH, CYPOR shRNA Plasmid (m): sc-35148-SH, CYPOR shRNA (h) Lentiviral Particles: sc-35147-V and CYPOR shRNA (m) Lentiviral Particles: sc-35148-V.


Molecular Weight of CYPOR: 76 kDa.

Positive Controls: MIA PaCa-2 cell lysate: sc-2285, Hep G2 cell lysate: sc-2227 or KNRK whole cell lysate: sc-2214.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 3) Immunohistochemistry: use ImmunoCruz™: sc-2053 or ABC: sc-2023 goat IgG Staining Systems.

DATA



CYPOR (C-19): sc-8091. Immunoperoxidase staining of formalin fixed, paraffin-embedded human liver tissue showing cytoplasmic staining of hepatocytes.

SELECT PRODUCT CITATIONS

- Ozkan, A., et al. 2004. Cytotoxicity of low dose epirubicin-HCl combined with lymphokine activated killer cells against hepatocellular carcinoma cell line hepatoma G₂. Turk. J. Med. Sci. 34: 11-19.
- Ozkan, A., et al. 2004. Epirubicin HCl toxicity in human-liver derived hepatoma G₂ cells. Pol. J. Pharmacol. 56: 435-444.
- Ozkan, A., et al. 2006. Protective effect of antioxidant enzymes against drug cytotoxicity in MCF7 cells. Exp. Oncol. 28: 86-88.
- Villaggio, B., et al. 2012. 1,25-dihydroxyvitamin D₃ downregulates aromatase expression and inflammatory cytokines in human macrophages. Clin. Exp. Rheumatol. 30: 934-938.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfation
Guaranteed

Try **CYPOR (F-10): sc-25270** or **CYPOR (F-2): sc-55477**, our highly recommended monoclonal alternatives to CYPOR (C-19).