

Factor X (MA-HFX): sc-81680

BACKGROUND

Hemostasis following tissue injury involves the deployment of essential plasma procoagulants (Prothrombin and Factors X, IX, V and VIII), which are involved in a blood coagulation cascade that leads to the formation of insoluble Fibrin clots and the promotion of platelet aggregation. Coagulation Factor X (Stuart Prower factor, FX, F10) is a vitamin K-dependent, single chain serine protease that is synthesized in the liver and circulates as an inactive precursor. The mature form of FX (FXA) is generated by Factor IX A- or Factor VII A-mediated cleavage at the tripeptide sequence, Arg-Lys-Arg, to yield a disulfide linked dimer. Together with the cofactor FVA and Ca²⁺ on the surface of platelets or endothelial cells, Factor X A coordinates as part of the prothrombinase complex, which mediates proteolysis of Prothrombin into active Thrombin. Mutations at the FX locus resulting in Factor X deficiencies can contribute to hemorrhagic diathesis.

REFERENCES

1. Davie, E.W. and Fujikawa, K. 1975. Basic mechanisms in blood coagulation. *Annu. Rev. Biochem.* 44: 799-829.
2. Di Scipio, R.G., Hermodson, M.A., Yates, S.G. and Davie, E.W. 1977. A comparison of human Prothrombin, Factor IX (Christmas factor), Factor X (Stuart factor), and Protein S. *Biochemistry* 16: 698-706.
3. Davie, E.W., Fujikawa, K. and Kisiel, W. 1991. The coagulation cascade: initiation, maintenance, and regulation. *Biochemistry* 30: 10363-10370.
4. Macedo-Ribeiro, S., Bode, W., Huber, R., Quinn-Allen, M.A., Kim, S.W., Ortel, T.L., Bourenkov, G.P., Bartunik, H.D., Stubbs, M.T., Kane, W.H. and Fuentes-Prior, P. 1999. Crystal structures of the membrane-binding C2 domain of human coagulation Factor V. *Nature* 402: 434-439.
5. Chambers, R.C., Leoni, P., Blanc-Brude, O.P., Wembridge, D.E. and Laurent, G.J. 2000. Thrombin is a potent inducer of connective tissue growth factor production via proteolytic activation of protease-activated receptor-1. *J. Biol. Chem.* 275: 35584-35591.
6. Yang, Y.H., Hwang, K.K., FitzGerald, J., Grossman, J.M., Taylor, M., Hahn, B.H. and Chen, P.P. 2006. Antibodies against the activated coagulation Factor X (FXA) in the antiphospholipid syndrome that interfere with the FXA inactivation by antithrombin. *J. Immunol.* 177: 8219-8225.
7. Todd, T., Perry, D.J., Hayman, E., Lawrence, K., Gattens, M. and Baglin, T. 2006. Severe Factor X deficiency due to a homozygous mutation (Cys364-Arg) that disrupts a disulphide bond in the catalytic domain. *Haemophilia* 12: 621-624.
8. Ndonwi, M., Broze, G.J., Jr., Agah, S., Schmidt, A.E. and Bajaj, S.P. 2007. Substitution of the Gla domain in Factor X with that of Protein C impairs its interaction with Factor VII A/tissue factor: lack of comparable effect by similar substitution in Factor IX. *J. Biol. Chem.* 282: 15632-15644.
9. Al-Hilali, A., Wulff, K., Abdel-Razeq, H., Saud, K.A., Al-Gaili, F. and Herrmann, F.H. 2007. Analysis of the novel Factor X gene mutation Glu51Lys in two families with Factor X-Riyadh anomaly. *Thromb. Haemost.* 97: 542-545.

CHROMOSOMAL LOCATION

Genetic locus: F10 (human) mapping to 13q34; F10 (mouse) mapping to 8 A1.1.

SOURCE

Factor X (MA-HFX) is a mouse monoclonal antibody raised against Factor X of human origin.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide, 0.1% gelatin and 5% glycerol.

APPLICATIONS

Sox-22 (86C2a) is recommended for detection of Sox-22 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) and immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)].

Suitable for use as control antibody for Sox-22 siRNA (h): sc-38435, Sox-22 shRNA Plasmid (h): sc-38435-SH and Sox-22 shRNA (h) Lentiviral Particles: sc-38435-V.

Molecular Weight of Factor X: 61 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-mouse IgG-HRP: sc-2005 (dilution range: 1:2000-1:32,000) or Cruz Marker™ compatible goat anti-mouse IgG-HRP: sc-2031 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.