

USP29 (C-20): sc-82083

BACKGROUND

The ubiquitin (Ub) pathway involves three sequential enzymatic steps that facilitate the conjugation of Ub and Ub-like molecules to specific protein substrates. Through the use of a wide range of enzymes that can add or remove ubiquitin, the Ub pathway controls many intracellular processes such as signal transduction, transcriptional activation and cell cycle progression. USP29 (ubiquitin specific peptidase 29), also known as HOM-TES-84/86, is a 922 amino acid deubiquitinating enzyme that participates in the Ub pathway. A member of the peptidase C19 family, the catalytic activity of USP29 involves a combination of the ubiquitin carboxyl-terminal thiolester and water to produce ubiquitin and a thiol. USP29 contains a cys box and a his box, which are characteristic of type-2 ubiquitin C-terminal hydrolases.

REFERENCES

1. D'Andrea, A. and Pellman, D. 1998. Deubiquitinating enzymes: a new class of biological regulators. *Crit. Rev. Biochem. Mol. Biol.* 33: 337-352.
2. Chung, C.H. and Baek, S.H. 1999. Deubiquitinating enzymes: their diversity and emerging roles. *Biochem. Biophys. Res. Commun.* 266: 633-640.
3. Kim, J., Noskov, V.N., Lu, X., Bergmann, A., Ren, X., Warth, T., Richardson, P., Kouprina, N. and Stubbs, L. 2000. Discovery of a novel, paternally expressed ubiquitin-specific processing protease gene through comparative analysis of an imprinted region of mouse chromosome 7 and human chromosome 19q13.4. *Genome Res.* 10: 1138-1147.
4. Kim, J., Bergmann, A., Wehri, E., Lu, X. and Stubbs, L. 2001. Imprinting and evolution of two Kruppel-type zinc-finger genes, ZIM3 and ZNF264, located in the PEG3/USP29 imprinted domain. *Genomics* 77: 91-98.
5. Online Mendelian Inheritance in Man, OMIM™. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 609546. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Kim, J., Bergmann, A., Choo, J.H. and Stubbs, L. 2007. Genomic organization and imprinting of the Peg3 domain in bovine. *Genomics* 90: 85-92.
7. Kim, J.D., Yu, S., Choo, J.H. and Kim, J. 2008. Two evolutionarily conserved sequence elements for Peg3/Usp29 transcription. *BMC Mol. Biol.* 9: 108.

CHROMOSOMAL LOCATION

Genetic locus: USP29 (human) mapping to 19q13.43.

SOURCE

USP29 (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of USP29 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-82083 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

USP29 (C-20) is recommended for detection of USP29 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other USP family members.

Suitable for use as control antibody for USP29 siRNA (h): sc-76833, USP29 shRNA Plasmid (h): sc-76833-SH and USP29 shRNA (h) Lentiviral Particles: sc-76833-V.

Molecular Weight of USP29: 104 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfation
Guaranteed

Try **USP29 (1A8): sc-517145**, our highly recommended monoclonal alternative to USP29 (C-20).