

MAML3 (C-19): sc-82220

BACKGROUND

MAML3 (mastermind-like protein 3) is a nuclear speckle protein that acts as a transcriptional coactivator for Notch receptors. The Notch signaling pathway influences cell fate by regulating the ability of precursor cells to properly respond to developmental signals. MAML3 is a member of the mastermind-like family of proteins that are human homologs of the *Drosophila melanogaster* mastermind protein. Through its N-terminal region, MAML3 interacts with the ankyrin repeats of the Notch proteins Notch 1, Notch 2, Notch 3 and Notch 4. This interaction leads to formation of a DNA-binding complex with the Notch proteins and RBP-J κ ; a complex that can then induce HES1 gene expression. While the N-terminal domain of MAML3 is essential for proper Notch binding, the C-terminal domain of MAML3 is essential for transcriptional activation. Due to its involvement in cell signaling and transcriptional activation, up-regulation of MAML3 is thought to be involved in oncogenesis.

REFERENCES

1. Wu, L., Sun, T., Kobayashi, K., Gao, P. and Griffin, J.D. 2002. Identification of a family of mastermind-like transcriptional coactivators for mammalian Notch receptors. *Mol. Cell. Biol.* 22: 7688-7700.
2. Lin, S.E., Oyama, T., Nagase, T., Harigaya, K. and Kitagawa, M. 2002. Identification of new human mastermind proteins defines a family that consists of positive regulators for Notch signaling. *J. Biol. Chem.* 277: 50612-50620.
3. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 608991. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Wu, L. and Griffin, J.D. 2004. Modulation of Notch signaling by mastermind-like (MAML) transcriptional co-activators and their involvement in tumorigenesis. *Semin. Cancer Biol.* 14: 348-356.
5. Katoh, M. and Katoh, M. 2006. Wnt antagonist, Dkk-2, is a Notch signaling target in intestinal stem cells: augmentation of a negative regulation system for canonical Wnt signaling pathway by the Notch-Dkk-2 signaling loop in primates. *Int. J. Mol. Med.* 19: 197-201.
6. William, D.A., Saitta, B., Gibson, J.D., Traas, J., Markov, V., Gonzalez, D.M., Sewell, W., Anderson, D.M., Pratt, S.C., Rappaport, E.F. and Kusumi, K. 2007. Identification of oscillatory genes in somitogenesis from functional genomic analysis of a human mesenchymal stem cell model. *Dev. Biol.* 305: 172-186.
7. Wu, L., Maillard, I., Nakamura, M., Pear, W.S. and Griffin, J.D. 2007. The transcriptional coactivator Maml1 is required for Notch2-mediated marginal zone B cell development. *Blood* 110: 3618-3623.

CHROMOSOMAL LOCATION

Genetic locus: MAML3 (human) mapping to 4q31.1; Maml3 (mouse) mapping to 3 C.

SOURCE

MAML3 (C-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of MAML3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-82220 X, 200 μ g/0.1 ml.

Blocking peptide available for competition studies, sc-82220 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MAML3 (C-19) is recommended for detection of MAML3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family member MAML2.

MAML3 (C-19) is also recommended for detection of MAML3 in additional species, including equine, porcine and avian.

Suitable for use as control antibody for MAML3 siRNA (h): sc-75744, MAML3 siRNA (m): sc-75745, MAML3 shRNA Plasmid (h): sc-75744-SH, MAML3 shRNA Plasmid (m): sc-75745-SH, MAML3 shRNA (h) Lentiviral Particles: sc-75744-V and MAML3 shRNA (m) Lentiviral Particles: sc-75745-V.

MAML3 (C-19) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of unprocessed MAML3 precursor: 150 kDa.

Molecular Weight of posttranslationally modified MAML3: 170 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

MONOS
Satisfaction
Guaranteed

Try **MAML3 (797C2a): sc-81104**, our highly recommended monoclonal alternative to MAML3 (C-19).