



# Pancreasin (M-20): sc-82276

## BACKGROUND

Serine proteases are important in many biological processes such as receptor activation, complement activation, coagulation and tissue remodeling. Pancreasin, also known as marapsin (MPN), channel activating protease 2-like protein (CAHP2) or protease, serine 27, is an N-glycosylated, secreted pancreatic tryptic serine peptidase and proteinase. Pancreasin is responsible for cleaving peptides after an Arginine residue and may play a role in regulating cell growth and migration. It can be inhibited by benzamidine and Leupeptin. Pancreasin is closely related to prostatic, Trypsin  $\gamma$ , Testisin and Trypsin  $\epsilon$ . These proteins share approximately 40% amino acid identity with trypsin  $\alpha$  and trypsin  $\beta$ . They contain Cysteine residues that may form a disulfide link between the propeptide and catalytic chain, a tryptic propeptide cleavage site and a C-terminal membrane anchor. Trypsin  $\epsilon$  and the human Pancreasin protein lack the characteristic C-terminal membrane anchor.

## REFERENCES

1. Bhagwandin, V.J., Hau, L.W., Mallen-St. Clair, J., Wolters, P.J. and Caughey, G.H. 2003. Structure and activity of human Pancreasin, a novel tryptic serine peptidase expressed primarily by the pancreas. *J. Biol. Chem.* 278: 3363-3371.
2. Tong, Z., Illek, B., Bhagwandin, V.J., Verghese, G.M. and Caughey, G.H. 2004. Prostatic, a membrane-anchored serine peptidase, regulates sodium currents in JME/CF15 cells, a cystic fibrosis airway epithelial cell line. *Am. J. Physiol. Lung Cell Mol. Physiol.* 287: L928-L935.
3. Verghese, G.M., Tong, Z.Y., Bhagwandin, V. and Caughey, G.H. 2004. Mouse Prostatic gene structure, promoter analysis, and restricted expression in lung and kidney. *Am. J. Respir. Cell Mol. Biol.* 30: 519-529.
4. Yasuda, S., Morokawa, N., Wong, G.W., Rossi, A., Madhusudhan, M.S., Sali, A., Askew, Y.S., Adachi, R., Silverman, G.A., Krilis, S.A. and Stevens, R.L. 2005. Urokinase-type plasminogen activator is a preferred substrate of the human epithelium serine protease Trypsin  $\epsilon$ /PRSS22. *Blood* 105: 3893-3901.
5. Wong, G.W. and Stevens, R.L. 2005. Identification of a subgroup of glycosylphosphatidylinositol-anchored trypsinases. *Biochem. Biophys. Res. Commun.* 336: 579-584.
6. Cal, S., Peinado, J.R., Llamazares, M., Quesada, V., Moncada-Pazos, A., Garabaya, C. and López-Otín, C. 2006. Identification and characterization of human Polyserine-3, a novel protein with tandem serine-protease domains in the same polypeptide chain. *BMC Biochem.* 7: 9-9.
7. Adachi, W., Ulanovsky, H., Li, Y., Norman, B., Davis, J. and Piatigorsky, J. 2006. Serial analysis of gene expression (SAGE) in the rat limbal and central corneal epithelium. *Invest. Ophthalmol. Vis. Sci.* 47: 3801-3810.
8. Stevens, R.L. and Adachi, R. 2007. Protease-proteoglycan complexes of mouse and human mast cells and importance of their trypsin  $\beta$ -heparin complexes in inflammation and innate immunity. *Immunol. Rev.* 217: 155-167.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## CHROMOSOMAL LOCATION

Genetic locus: Prss27 (mouse) mapping to 17 A3.3.

## SOURCE

Pancreasin (M-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Pancreasin of mouse origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-82276 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

Pancreasin (M-20) is recommended for detection of Pancreasin of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Pancreasin siRNA (m): sc-72313, Pancreasin shRNA Plasmid (m): sc-72313-SH and Pancreasin shRNA (m) Lentiviral Particles: sc-72313-V.

Molecular Weight of Pancreasin: 32 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotting A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.