

Manic Fringe (N-19): sc-8238

BACKGROUND

Three mammalian fringe family members, Manic, Radical and Lunatic Fringe, have been identified as proteins related to *Drosophila* fringe, a protein involved in development. Fringe proteins act upstream of the Notch signaling pathway and are involved in boundary determination during segmentation. Each mammalian fringe displays different patterns of expression, though all are expressed in the mouse embryo as well as in many adult tissues. Radical fringe plays a key role in the development of the limb bud. Lunatic fringe is required for normal somite segmentation and patterning and is thought to be a target of the molecular clock. Manic fringe, also involved in somatic development, has been shown to render mouse NIH/3T3 cells tumorigenic in SCID mice.

REFERENCES

1. May, W.A., et al. 1997. EWS/FLI1-induced manic fringe renders NIH/3T3 cells tumorigenic. *Nat. Genet.* 17: 495-497.
2. Laufer, E., et al. 1997. Expression of Radical fringe in limb-bud ectoderm regulates apical ectodermal ridge formation. *Nature* 386: 366-373.
3. Johnston, S.H., et al. 1997. A family of mammalian Fringe genes implicated in boundary determination and the Notch pathway. *Development* 124: 2245-2254.
4. May, W.A., et al. 1997. EWS/FLI1-induced manic fringe renders NIH 3T3 cells tumorigenic. *Nat. Genet.* 17: 495-497.
5. Cohen, B., et al. 1997. Fringe boundaries coincide with Notch-dependent patterning centres in mammals and alter Notch-dependent development in *Drosophila*. *Nat. Genet.* 16: 283-288.
6. Thelu, J., et al. 1998. Differential expression pattern of the three fringe genes is associated with epidermal differentiation. *J. Invest. Dermatol.* 111: 903-906.
7. Evrard, Y.A., et al. 1998. Lunatic fringe is an essential mediator of somite segmentation and patterning. *Nature* 394: 377-381.
8. McGrew, M.J., et al. 1998. The lunatic fringe gene is a target of the molecular clock linked to somite segmentation in avian embryos. *Curr. Biol.* 8: 979-982.

CHROMOSOMAL LOCATION

Genetic locus: MFNG (human) mapping to 22q13.1.

SOURCE

Manic Fringe (N-19) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of Manic Fringe of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-8238 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Manic Fringe (N-19) is recommended for detection of Manic Fringe of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Manic Fringe siRNA (h): sc-39492, Manic Fringe shRNA Plasmid (h): sc-39492-SH and Manic Fringe shRNA (h) Lentiviral Particles: sc-39492-V.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

SELECT PRODUCT CITATIONS

1. Chakrabarti, O., et al. 2004. Human papillomavirus type 16 E6 amino acid 83 variants enhance E6-mediated MAPK signaling and differentially regulate tumorigenesis by notch signaling and oncogenic Ras. *J. Virol.* 78: 5934-5945.
2. Veeraraghavalu, K. 2004. Papillomavirus-mediated neoplastic progression is associated with reciprocal changes in JAGGED1 and Manic Fringe expression linked to notch activation. *J. Virol.* 78: 8687-8700.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
 Satisfaction
 Guaranteed

Try **Manic Fringe (2B11): sc-517162**, our highly recommended monoclonal alternative to Manic Fringe (N-19).