



Zimp10 (C-20): sc-82437

BACKGROUND

Zimp10, also known as ZMIZ1 (zinc finger, MIZ-type containing 1), MIZ or RAI17, is a 1,067 amino acid protein that localizes to both the nucleus and the cytoplasm and contains one SP-RING-type zinc finger. Expressed in ovary, prostate, testis and spleen, Zimp10 interacts with AR (androgen receptor) and functions to increase the ligand-dependent transcriptional activity of AR and promote the sumoylation of AR, an event which is necessary for AR stimulation and overall activity. Multiple isoforms of Zimp10 exist due to alternative splicing events. The gene encoding Zimp10 maps to human chromosome 10, which houses over 1,200 genes and comprises nearly 4.5% of the human genome. Defects in some of the genes that map to chromosome 10 are associated with Charcot-Marie tooth disease, Jackson-Weiss syndrome, Usher syndrome, nonsyndromic deafness, Wolman's syndrome, Cowden syndrome, multiple endocrine neoplasia type 2 and porphyria.

REFERENCES

- Nagase, T., Ishikawa, K., Kikuno, R., Hirosawa, M., Nomura, N. and Ohara, O. 1999. Prediction of the coding sequences of unidentified human genes. XV. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. DNA Res. 6: 337-345.
- Sharma, M., Li, X., Wang, Y., Zarnegar, M., Huang, C.Y., Palvimo, J.J., Lim, B. and Sun, Z. 2003. hZimp10 is an androgen receptor co-activator and forms a complex with SUMO-1 at replication foci. EMBO J. 22: 6101-6114.
- Li, X., Thyssen, G., Beliakoff, J. and Sun, Z. 2006. The novel PIAS-like protein hZimp10 enhances Smad transcriptional activity. J. Biol. Chem. 281: 23748-23756.
- Lee, J., Beliakoff, J. and Sun, Z. 2007. The novel PIAS-like protein hZimp10 is a transcriptional co-activator of the p53 tumor suppressor. Nucleic Acids Res. 35: 4523-4534.
- Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 607159. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Soler, G., Radford-Weiss, I., Ben-Abdelali, R., Mahlaoui, N., Ponceau, J.F., Macintyre, E.A., Vekemans, M., Bernard, O.A. and Romana, S.P. 2008. Fusion of ZMIZ1 to ABL1 in a B cell acute lymphoblastic leukaemia with a t(9;10)(q34;q22.3) translocation. Leukemia 22: 1278-1280.
- Beliakoff, J., Lee, J., Ueno, H., Aiyer, A., Weissman, I.L., Barsh, G.S., Cardiff, R.D. and Sun, Z. 2008. The PIAS-like protein Zimp10 is essential for embryonic viability and proper vascular development. Mol. Cell. Biol. 28: 282-292.

CHROMOSOMAL LOCATION

Genetic locus: ZMIZ1 (human) mapping to 10q22.3; Zmiz1 (mouse) mapping to 14 A3.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

SOURCE

Zimp10 (C-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Zimp10 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-82437 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Zimp10 (C-20) is recommended for detection of Zimp10 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Zimp10 siRNA (h): sc-76960, Zimp10 siRNA (m): sc-76961, Zimp10 shRNA Plasmid (h): sc-76960-SH, Zimp10 shRNA Plasmid (m): sc-76961-SH, Zimp10 shRNA (h) Lentiviral Particles: sc-76960-V and Zimp10 shRNA (m) Lentiviral Particles: sc-76961-V.

Molecular Weight of Zimp10: 130 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.