

HGSNAT (K-15): sc-82623

BACKGROUND

HGSNAT (heparan- α -glucosaminide N-acetyltransferase), also known as HGNAT, MPS3C or TMEM76, is a 663 amino acid multi-pass membrane protein that localizes to the lysosome and is expressed in a wide variety of tissues, including liver, lung, placenta, heart and skeletal muscle. Functioning as a lysosomal acetyltransferase, HGSNAT acetylates the non-reducing terminal α -glucosamine residue of intralysosomal heparin or heparan sulfate, effectively producing a substrate for luminal α -N-acetyl glucosaminidase. Defects in the gene encoding HGSNAT are the cause of mucopolysaccharidosis type 3C (MPS3C), an autosomal recessive lysosomal storage disease that is caused by impaired heparan sulfate degradation and is characterized by degeneration of the central nervous system.

REFERENCES

1. Bame, K.J. and Rome, L.H. 1986. Genetic evidence for transmembrane acetylation by lysosomes. *Science* 233: 1087-1089.
2. Meikle, P.J., Whittle, A.M. and Hopwood, J.J. 1995. Human acetyl-coenzyme A: α -glucosaminide N-acetyltransferase. Kinetic characterization and mechanistic interpretation. *Biochem. J.* 308: 327-333.
3. Ausseil, J., Loredó-Ostí, J.C., Verner, A., Darmond-Zwaig, C., Maire, I., Poorthuis, B., van Diggelen, O.P., Hudson, T.J., Fujiwara, T.M., Morgan, K. and Pshezhetsky, A.V. 2004. Localisation of a gene for mucopolysaccharidosis IIIC to the pericentromeric region of chromosome 8. *J. Med. Genet.* 41: 941-945.
4. Fan, X., Zhang, H., Zhang, S., Bagshaw, R.D., Tropak, M.B., Callahan, J.W. and Mahuran, D.J. 2006. Identification of the gene encoding the enzyme deficient in mucopolysaccharidosis IIIC (Sanfilippo disease type C). *Am. J. Hum. Genet.* 79: 738-744.
5. Hřebíček, M., Mrázová, L., Seyrantepe, V., Durand, S., Roslin, N.M., Nosková, L., Hartmannová, H., Ivánek, R., Cízková, A., Poupětová, H., Sikora, J., Urinová, J., Stranecký, V., Zeman, J., Lepage, P., Roquis, D., et al. 2006. Mutations in TMEM76* cause mucopolysaccharidosis IIIC (Sanfilippo C syndrome). *Am. J. Hum. Genet.* 79: 807-819.
6. Fedele, A.O., Filocamo, M., Di Rocco, M., Sersale, G., Lübke, T., di Natale, P., Cosma, M.P. and Ballabio, A. 2007. Mutational analysis of the HGSNAT gene in Italian patients with mucopolysaccharidosis IIIC (Sanfilippo C syndrome). Mutation in brief #959. Online. *Hum. Mutat.* 28: 523.
7. Coutinho, M.F., Lacerda, L., Prata, M.J., Ribeiro, H., Lopes, L., Ferreira, C. and Alves, S. 2008. Molecular characterization of Portuguese patients with mucopolysaccharidosis IIIC: two novel mutations in the HGSNAT gene. *Clin. Genet.* 74: 194-195.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 610453. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: HGSNAT (human) mapping to 8p11.21; Hgsnat (mouse) mapping to 8 A2.

SOURCE

HGSNAT (K-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of HGSNAT of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-82623 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HGSNAT (K-15) is recommended for detection of HGSNAT of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

HGSNAT (K-15) is also recommended for detection of HGSNAT in additional species, including equine, canine, porcine and avian.

Suitable for use as control antibody for HGSNAT siRNA (h): sc-75251, HGSNAT siRNA (m): sc-75252, HGSNAT shRNA Plasmid (h): sc-75251-SH, HGSNAT shRNA Plasmid (m): sc-75252-SH, HGSNAT shRNA (h) Lentiviral Particles: sc-75251-V and HGSNAT shRNA (m) Lentiviral Particles: sc-75252-V.

Molecular Weight of HGSNAT: 73 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.