

HoxC9 (G-14): sc-82914

BACKGROUND

The Hox proteins are a family of transcription factors that play a role in development and cellular differentiation by regulating downstream target genes. Specifically, the Hox proteins direct DNA-protein and protein-protein interactions that assist in determining the morphologic features associated with the anterior-posterior body axis. Hox proteins are involved in controlling axial patterning, leukemias and hereditary malformations. HoxC9 (homeobox protein HoxC9), also known as HOX3 or HOX3B, is a member of the Abd-B homeobox (Hox) family. It is a 260 amino acid long nuclear protein that contains one homeobox DNA-binding domain. HoxC9 plays a role in the regulation of development, providing cells with positional identities on the anterior-posterior body axis. In addition, HoxC9 is expressed in esophageal cancer cells and may be involved in cancer development.

REFERENCES

1. Redline, R.W., Hudock, P., MacFee, M. and Patterson, P. 1994. Expression of Abd-B-type homeobox genes in human tumors. *Lab. Invest.* 71: 663-670.
2. Miano, J.M., Firulli, A.B., Olson, E.N., Hara, P., Giachelli, C.M. and Schwartz, S.M. 1996. Restricted expression of homeobox genes distinguishes fetal from adult human smooth muscle cells. *Proc. Natl. Acad. Sci. USA* 93: 900-905.
3. Ponsuksili, S., Wimmers, K., Adjaye, J. and Schellander, K. 2001. Expression of homeobox-containing genes in cDNA libraries derived from cattle oocytes and preimplantation stage embryo. *Mol. Reprod. Dev.* 60: 297-301.
4. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 142970. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. López, R., Garrido, E., Vázquez, G., Piña, P., Pérez, C., Alvarado, I. and Salcedo, M. 2006. A subgroup of Hox Abd-B gene is differentially expressed in cervical cancer. *Int. J. Gynecol. Cancer* 16: 1289-1296.
6. Huntriss, J., Hinkins, M. and Picton, H.M. 2006. cDNA cloning and expression of the human NOBOX gene in oocytes and ovarian follicles. *Mol. Hum. Reprod.* 12: 283-289.
7. Okamoto, O.K., Oba-Shinjo, S.M., Lopes, L. and Nagahashi Marie, S.K. 2007. Expression of HOXC9 and E2F2 are upregulated in CD133+ cells isolated from human astrocytomas and associate with transformation of human astrocytes. *Biochim. Biophys. Acta* 1769: 437-442.

CHROMOSOMAL LOCATION

Genetic locus: HOXC9 (human) mapping to 12q13.13; Hoxc9 (mouse) mapping to 15 F3.

SOURCE

HoxC9 (G-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of HoxC9 of human origin.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-82914 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

HoxC9 (G-14) is recommended for detection of HoxC9 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other Hox family members.

HoxC9 (G-14) is also recommended for detection of HoxC9 in additional species, including bovine.

Suitable for use as control antibody for HoxC9 siRNA (h): sc-75289, HoxC9 siRNA (m): sc-75290, HoxC9 shRNA Plasmid (h): sc-75289-SH, HoxC9 shRNA Plasmid (m): sc-75290-SH, HoxC9 shRNA (h) Lentiviral Particles: sc-75289-V and HoxC9 shRNA (m) Lentiviral Particles: sc-75290-V.

HoxC9 (G-14) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.

Molecular Weight of HoxC9: 29 kDa.

Positive Controls: HT-1080 whole cell lysate: sc-364183.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **HoxC9 (HOXCA6E6): sc-81100**, our highly recommended monoclonal alternative to HoxC9 (G-14).