

STCH (C-20): sc-83335

BACKGROUND

The heat shock proteins (HSPs) comprise a group of highly conserved, abundantly expressed proteins with diverse functions, which include the assembly and sequestering of multiprotein complexes, transportation of nascent polypeptide chains across cellular membranes and regulation of protein folding. Heat shock proteins (also known as molecular chaperones) fall into six general families: HSP 90, HSP 70, HSP 60, the low molecular weight HSPs, the immunophilins and the HSP 105 family. STCH, also known as HSPA13 (heat shock 70 kDa protein, member 13) or microsomal stress 70 protein ATPase core, is a 471 amino acid protein belonging to the HSP 70 protein family. Localized to the endoplasmic reticulum, STCH has peptide-independent ATPase activity. Fundamentally expressed in all tissues, STCH has been found to interact with PLIC-1 and PLIC-2, proteins involved in the signaling connection between the membrane receptors for Thrombospondin and the cytoskeleton.

REFERENCES

- Georgopoulos, C. and Welch, W.J. 1993. Role of the major heat shock proteins as molecular chaperones. *Annu. Rev. Cell Biol.* 9: 601-634.
- Otterson, G.A., et al. 1994. Stch encodes the "TPase core" of a microsomal stress 70 protein. *EMBO J.* 13: 1216-1225.
- Todd, M.J., et al. 1994. Dynamics of the chaperonin ATPase cycle: implications for facilitated protein folding. *Science* 265: 659-666.
- Brodsky, G., et al. 1995. Localization of STCH to human chromosome 21q11.1. *Genomics* 30: 627-628.
- Online Mendelian Inheritance in Man, OMIM[™]. 1999. Johns Hopkins University, Baltimore, MD. MIM Number: 601100. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Kaye, F.J., et al. 2000. A family of ubiquitin-like proteins binds the ATPase domain of HSP 70-like Stch. *FEBS Lett.* 467: 348-355.
- Chouchane, L., et al. 2001. Genetic variation in the stress protein HSP 70-2 gene is highly associated with obesity. *Int. J. Obes. Relat. Metab. Disord.* 25: 462-466.
- Aoki, M., et al. 2005. A genetic variant in the gene encoding the stress 70 protein chaperone family member STCH is associated with gastric cancer in the Japanese population. *Biochem. Biophys. Res. Commun.* 335: 566-574.
- Yamagata, N., et al. 2008. Stomach cancer-derived del223V-226L mutation of the STCH gene causes loss of sensitization to TRAIL-mediated apoptosis. *Biochem. Biophys. Res. Commun.* 376: 499-503.

CHROMOSOMAL LOCATION

Genetic locus: HSPA13 (human) mapping to 21q11.2; Hspa13 (mouse) mapping to 16 C3.1.

SOURCE

STCH (C-20) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping at the C-terminus of STCH of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-83335 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

STCH (C-20) is recommended for detection of STCH of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

STCH (C-20) is also recommended for detection of STCH in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for STCH siRNA (h): sc-91484, STCH siRNA (m): sc-153884, STCH shRNA Plasmid (h): sc-91484-SH, STCH shRNA Plasmid (m): sc-153884-SH, STCH shRNA (h) Lentiviral Particles: sc-91484-V and STCH shRNA (m) Lentiviral Particles: sc-153884-V.

Molecular Weight of STCH: 60 kDa.

Positive Controls: U-937 cell lysate: sc-2239, Hep G2 cell lysate: sc-2227 or K-562 whole cell lysate: sc-2203.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **STCH (A-11): sc-398297**, our highly recommended monoclonal alternative to STCH (C-20).