

ZnT-9 (E-17): sc-83505

BACKGROUND

Zinc, an essential element required for cell proliferation and differentiation, plays a role in a diverse array of cellular functions, including acting as a cofactor for numerous enzymes and transcription factors and as a neuroregulator. The zinc transporter (ZnT) or solute carrier 30 (SLC30) family regulates the supply of zinc within cells, and its members are characterized by containing six membrane-spanning domains, a large histidine-rich intracellular loop, and a C-terminal tail. ZnT proteins also belong to the cation diffusion facilitator (CDF) transporter family of metal ion transporters. ZnT-9, also known as HUEL (human embryonic lung protein), GAC63 (GRIP1-associated co-activator 1) or SLC30 member 9, displays ubiquitous expression in fetal and adult tissues as well as cancer cell lines. ZnT-9 localizes to the cytoplasm and is translocated to the nucleus during S phase. ZnT-9 has the lowest homology with the other zinc transporters and may function as a DNA-binding protein.

REFERENCES

- Palmiter, R.D. and Huang, L. 2004. Efflux and compartmentalization of zinc by members of the SLC30 family of solute carriers. *Pflugers Arch.* 447: 744-751.
- Seve, M., Chimienti, F., Devergnas, S. and Favier, A. 2004. In silico identification and expression of SLC30 family genes: an expressed sequence tag data mining strategy for the characterization of zinc transporters' tissue expression. *BMC Genomics* 5: 32-32.
- González-Guerrero, M., Azcón-Aguilar, C., Mooney, M., Valderas, A., MacDiarmid, C.W., Eide, D.J. and Ferrol, N. 2005. Characterization of a *Glomus* intraradices gene encoding a putative Zn transporter of the cation diffusion facilitator family. *Fungal Genet. Biol.* 42: 130-140.
- Chimienti, F., Favier, A. and Seve, M. 2005. ZnT-8, a pancreatic β cell-specific zinc transporter. *Biometals* 18: 313-317.
- Kambe, T., Suzuki, T., Nagao, M. and Yamaguchi-Iwai, Y. 2006. Sequence similarity and functional relationship among eukaryotic ZIP and CDF transporters. *Genomics Proteomics Bioinformatics* 4: 1-9.
- Falcón-Perez, J.M. and Dell'Angelica, E.C. 2007. Zinc transporter 2 (SLC30A2) can suppress the vesicular zinc defect of adaptor protein 3-depleted fibroblasts by promoting zinc accumulation in lysosomes. *Exp. Cell Res.* 313: 1473-1483.
- Smidt, K., Pedersen, S.B., Brock, B., Schmitz, O., Fisker, S., Bendix, J., Wogensen, L. and Rungby, J. 2007. Zinc-transporter genes in human visceral and subcutaneous adipocytes: lean versus obese. *Mol. Cell. Endocrinol.* 264: 68-73.
- Montanini, B., Blaudez, D., Jeandroz, S., Sanders, D. and Chalot, M. 2007. Phylogenetic and functional analysis of the Cation Diffusion Facilitator (CDF) family: improved signature and prediction of substrate specificity. *BMC Genomics* 8: 107-107.

CHROMOSOMAL LOCATION

Genetic locus: SLC30A9 (human) mapping to 4p13; Slc30a9 (mouse) mapping to 5 C3.1.

SOURCE

ZnT-9 (E-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of ZnT-9 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-83505 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ZnT-9 (E-17) is recommended for detection of ZnT-9 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other ZnT family members.

ZnT-9 (E-17) is also recommended for detection of ZnT-9 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for ZnT-9 siRNA (h): sc-77015, ZnT-9 siRNA (m): sc-77016, ZnT-9 shRNA Plasmid (h): sc-77015-SH, ZnT-9 shRNA Plasmid (m): sc-77016-SH, ZnT-9 shRNA (h) Lentiviral Particles: sc-77015-V and ZnT-9 shRNA (m) Lentiviral Particles: sc-77016-V.

Molecular Weight of ZnT-9: 64 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.