

3 β -HSD7 (A-13): sc-83772

BACKGROUND

3 β -hydroxysteroid dehydrogenase (3 β -HSD), also known as HSD3B1 or HSDB3, is a bifunctional enzyme that plays a crucial role in the synthesis of all classes of hormonal steroids. 3 β -hydroxysteroid dehydrogenase type 7 (3 β -HSD7) is a 369 amino acid protein belonging to the 3 β -HSD family. Localized to the membrane of the endoplasmic reticulum, 3 β -HSD7 plays an active role in the initial stages of bile synthesis from cholesterol. 3 β -HSD7 catalyzes the oxidative conversion of δ 5-ene-3 β -hydroxy steroid and is active against 7 α hydroxylated sterol substrates. Mutations in the gene encoding 3 β -HSD7 are associated with congenital bile acid synthesis defect, which leads to neonatal cholestasis.

REFERENCES

- Cheng, J.B., Jacquemin, E., Gerhardt, M., Nazer, H., Cresteil, D., Heubi, J.E., Setchell, K.D. and Russell, D.W. 2003. Molecular genetics of 3 β -hydroxy- δ 5-C27-steroid oxidoreductase deficiency in 16 patients with loss of bile acid synthesis and liver disease. *J. Clin. Endocrinol. Metab.* 88: 1833-1841.
- Rahier, A., Darnet, S., Bouvier, F., Camara, B. and Bard, M. 2006. Molecular and enzymatic characterizations of novel bifunctional 3 β -hydroxysteroid dehydrogenases/C4 decarboxylases from *Arabidopsis thaliana*. *J. Biol. Chem.* 281: 27264-27277.
- Online Mendelian Inheritance in Man, OMIM[™]. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 607764. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Fischler, B., Bodin, K., Stjernman, H., Olin, M., Hansson, M., Sjövall, J. and Björkhem, I. 2007. Cholestatic liver disease in adults may be due to an inherited defect in bile acid biosynthesis. *J. Intern. Med.* 262: 254-262.
- Shea, H.C., Head, D.D., Setchell, K.D. and Russell, D.W. 2007. Analysis of HSD3B7 knockout mice reveals that a 3 α -hydroxyl stereochemistry is required for bile acid function. *Proc. Natl. Acad. Sci. USA* 104: 11526-11533.
- Claahsen-van der Grinten, H.L., Duthoi, K., Otten, B.J., d'Ancona, F.C., Hulsbergen-vd Kaa, C.A. and Hermus, A.R. 2008. An adrenal rest tumour in the perirenal region in a patient with congenital adrenal hyperplasia due to congenital 3 β -hydroxysteroid dehydrogenase deficiency. *Eur. J. Endocrinol.* 159: 489-491.
- Thomas, J.L., Mack, V.L., Glow, J.A., Moshkelani, D., Terrell, J.R. and Bucholtz, K.M. 2008. Structure/function of the inhibition of human 3 β -hydroxysteroid dehydrogenase type 1 and type 2 by trilostane. *J. Steroid Biochem. Mol. Biol.* 111: 66-73.
- Tripodi, G., Citterio, L., Kouznetsova, T., Lanzani, C., Florio, M., Modica, R., Messaggio, E., Hamlyn, J.M., Zagato, L., Bianchi, G., Staessen, J.A. and Manunta, P. 2009. Steroid biosynthesis and renal excretion in human essential hypertension: association with blood pressure and endogenous ouabain. *Am. J. Hypertens.* 22: 357-363.
- Rahier, A., Bergdoll, M., Genot, G., Bouvier, F. and Camara, B. 2009. Homology modeling and site-directed mutagenesis reveal catalytic key amino acids of 3 β -hydroxysteroid-dehydrogenase/C4-decarboxylase from *Arabidopsis thaliana*. *Plant Physiol.* 149: 1872-1886.

CHROMOSOMAL LOCATION

Genetic locus: HSD3B7 (human) mapping to 16p11.2; Hsd3b7 (mouse) mapping to 7 F3.

SOURCE

3 β -HSD7 (A-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of 3 β -HSD7 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-83772 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

3 β -HSD7 (A-13) is recommended for detection of 3 β -HSD7 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with 3 β -HSD.

3 β -HSD7 (A-13) is also recommended for detection of 3 β -HSD7 in additional species, including canine, bovine and porcine.

Suitable for use as control antibody for 3 β -HSD7 siRNA (h): sc-72402, 3 β -HSD7 siRNA (m): sc-72403, 3 β -HSD7 shRNA Plasmid (h): sc-72402-SH, 3 β -HSD7 shRNA Plasmid (m): sc-72403-SH, 3 β -HSD7 shRNA (h) Lentiviral Particles: sc-72402-V and 3 β -HSD7 shRNA (m) Lentiviral Particles: sc-72403-V.

Molecular Weight of 3 β -HSD7: 41 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.