

EBP (N-15): sc-83791

BACKGROUND

EBP (emopamil-binding protein), also known as CPX, CHO2, CPXD or CDPX2, is a 230 amino acid multi-pass membrane protein that localizes to the endoplasmic reticulum (ER) and is involved in steroid biosynthesis. Functioning to catalyze the conversion of Δ^8 -sterols to their corresponding Δ^7 -isomers, EBP plays an important role in drug transport and cholesterol metabolism within the ER. Defects in the gene encoding EBP are the cause of chondrodysplasia punctata X-linked dominant type 2 (CDPX2), a heterogeneous disorder that is caused by defective cholesterol biosynthesis. CDPX2 is characterized by punctiform calcification of the bones leading to linear ichthyosis, cataracts and short stature.

REFERENCES

- Hanner, M., Moebius, F.F., Weber, F., Grabner, M., Striessnig, J. and Glossmann, H. 1995. Phenylalkylamine Ca^{2+} antagonist binding protein. Molecular cloning, tissue distribution, and heterologous expression. *J. Biol. Chem.* 270: 7551-7557.
- Derry, J.M., Gormally, E., Means, G.D., Zhao, W., Meindl, A., Kelley, R.I., Boyd, Y. and Herman, G.E. 1999. Mutations in a Δ^8 - Δ^7 sterol isomerase in the tattered mouse and X-linked dominant chondrodysplasia punctata. *Nat. Genet.* 22: 286-290.
- Braverman, N., Lin, P., Moebius, F.F., Obie, C., Moser, A., Glossmann, H., Wilcox, W.R., Rimoin, D.L., Smith, M., Kratz, L., Kelley, R.I. and Valle, D. 1999. Mutations in the gene encoding 3β -hydroxysteroid- Δ^8 , Δ^7 -isomerase cause X-linked dominant Conradi-Hünermann syndrome. *Nat. Genet.* 22: 291-294.
- Has, C., Bruckner-Tuderman, L., Müller, D., Floeth, M., Folkers, E., Donnai, D. and Traupe, H. 2000. The Conradi-Hünermann-Happle syndrome (CDPX2) and emopamil binding protein: novel mutations, and somatic and gonadal mosaicism. *Hum. Mol. Genet.* 9: 1951-1955.
- Becker, K., Csikós, M., Horváth, A. and Kárpáti, S. 2001. Identification of a novel mutation in 3β -hydroxysteroid- Δ^8 - Δ^7 -isomerase in a case of Conradi-Hünermann-Happle syndrome. *Exp. Dermatol.* 10: 286-289.
- Moebius, F.F., Fitzky, B.U., Wietzorrek, G., Haidekker, A., Eder, A. and Glossmann, H. 2003. Cloning of an emopamil-binding protein (EBP)-like protein that lacks sterol Δ^8 - Δ^7 isomerase activity. *Biochem. J.* 374: 229-237.

CHROMOSOMAL LOCATION

Genetic locus: EBP (human) mapping to Xp11.23.

SOURCE

EBP (N-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the N-terminus of EBP of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-83791 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

EBP (N-15) is recommended for detection of EBP of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

EBP (N-15) is also recommended for detection of EBP in additional species, including canine.

Suitable for use as control antibody for EBP siRNA (h): sc-77218, EBP shRNA Plasmid (h): sc-77218-SH and EBP shRNA (h) Lentiviral Particles: sc-77218-V.


Molecular Weight of EBP: 26 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227, Jurkat whole cell lysate: sc-2204 or HeLa whole cell lysate: sc-2200.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



EBP (N-15): sc-83791. Western blot analysis of EBP expression in Jurkat (A), HeLa (B) and Hep G2 (C) whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.