

# AKAP 220 (G-20): sc-83965

## BACKGROUND

The type II cAMP-protein kinase (PKA) is a multifunctional kinase with a broad range of substrates. Specificity of PKA signaling is thought to be mediated by the compartmentalization of the kinase to specific sites within the cell. To maintain this specific localization, the R subunit (RII) of PKA interacts with specific RII-anchoring proteins. This family of proteins has been designated A-kinase anchoring proteins (AKAP). Members of this family, including MAP2 (microtubule-associated protein 2), neuronally expressed AKAP 79 and AKAP 150, and the DNA binding AKAP 95, display differential tissue specificity and localization. AKAP 220, may play a role in cAMP-responsive peroxisomal events by targeting type II PKA.

## REFERENCES

1. Scott, J.D., Stofko, R.E., McDonald, J.R., Comer, J.D., Vitalis, E.A. and Mangili, J.A. 1990. Type II regulatory subunit dimerization determines the subcellular localization of the cAMP-dependent protein kinase. *J. Biol. Chem.* 265: 21561-21566.
2. Carr, D.W., Stofko-Hahn, R.E., Fraser, I.D., Cone, R.D. and Scott, J.D. 1992. Localization of the cAMP-dependent protein kinase to the postsynaptic densities by A-kinase anchoring proteins. Characterization of AKAP 79. *J. Biol. Chem.* 267: 16816-16823.
3. Coghlan, V.M., Bergeson, S.E., Langeberg, L., Nilaver, G. and Scott, J.D. 1993. A-kinase anchoring proteins: a key to selective activation of cAMP-responsive events? *Mol. Cell. Biochem.* 127: 309-319.
4. Coghlan, V.M., Langeberg, L.K., Fernandez, A., Lamb, N.J. and Scott, J.D. 1994. Cloning and characterization of AKAP 95, a nuclear protein that associates with the regulatory subunit of type II cAMP-dependent protein kinase. *J. Biol. Chem.* 269: 7658-7665.
5. Coghlan, V.M., Perrino, B.A., Howard, M., Langeberg, L.K., Hicks, J.B., Gallatin, W.M. and Scott, J.D. 1995. Association of protein kinase A and protein phosphatase 2B with a common anchoring protein. *Science* 267: 108-111.
6. Lester, L.B., Coghlan, V.M., Nauert, B. and Scott, J.D. 1996. Cloning and characterization of a novel A-kinase anchoring protein. AKAP 220, association with testicular peroxisomes. *J. Biol. Chem.* 271: 9460-9465.

## CHROMOSOMAL LOCATION

Genetic locus: AKAP11 (human) mapping to 13q14.11; Akap11 (mouse) mapping to 14 D3.

## SOURCE

AKAP 220 (G-20) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping near the C-terminus of AKAP 220 of human origin.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-83965 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

AKAP 220 (G-20) is recommended for detection of AKAP 220 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

AKAP 220 (G-20) is also recommended for detection of AKAP 220 in additional species, including equine and bovine.

Suitable for use as control antibody for AKAP 220 siRNA (h): sc-105049, AKAP 220 siRNA (m): sc-40304, AKAP 220 shRNA Plasmid (h): sc-105049-SH, AKAP 220 shRNA Plasmid (m): sc-40304-SH, AKAP 220 shRNA (h) Lentiviral Particles: sc-105049-V and AKAP 220 shRNA (m) Lentiviral Particles: sc-40304-V.

Molecular Weight of AKAP 220: 220 kDa.

Positive Controls: Rat brain extract: sc-2392, rat testis extract: sc-2400 or KNRK whole cell lysate: sc-2214.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.