

MIPEP (S-18): sc-84203

BACKGROUND

MIPEP (mitochondrial intermediate peptidase) is a 713 amino acid protein that performs the final step in processing nuclear-encoded proteins that are targeted to the mitochondrial inner membrane or matrix. MIPEP is primarily involved in the maturation of oxidative phosphorylation-related proteins and is accordingly expressed at highest levels in tissues that have high oxygen consumption, such as skeletal muscle, brain, heart and pancreas. Significantly, human MIPEP is the functional homolog of yeast YMIP, a protein that promotes mitochondrial iron uptake and is regulated by the yeast homolog of Frataxin, which is required for mitochondrial iron efflux. Due to the likely regulatory effect that human Frataxin has on it, MIPEP may be implicated in the clinical manifestations of Friedreich ataxia, an autosomal recessive neurodegenerative disease caused by deficiency of Frataxin.

REFERENCES

1. Isaya, G., Kalousek, F. and Rosenberg, L.E. 1992. Sequence analysis of rat mitochondrial intermediate peptidase: similarity to zinc metallopeptidases and to a putative yeast homologue. *Proc. Natl. Acad. Sci. USA* 89: 8317-8321.
2. Chew, A., Buck, E.A., Peretz, S., Sirugo, G., Rinaldo, P. and Isaya, G. 1997. Cloning, expression, and chromosomal assignment of the human mitochondrial intermediate peptidase gene (MIPEP). *Genomics* 40: 493-496.
3. Branda, S.S., Yang, Z.Y., Chew, A. and Isaya, G. 1999. Mitochondrial intermediate peptidase and the yeast Frataxin homolog together maintain mitochondrial iron homeostasis in *Saccharomyces cerevisiae*. *Hum. Mol. Genet.* 8: 1099-1110.
4. Chew, A., Sirugo, G., Alsobrook, J.P. and Isaya, G. 2000. Functional and genomic analysis of the human mitochondrial intermediate peptidase, a putative protein partner of Frataxin. *Genomics* 65: 104-112.
5. Cavadini, P., Gellera, C., Patel, P.I. and Isaya, G. 2000. Human Frataxin maintains mitochondrial iron homeostasis in *Saccharomyces cerevisiae*. *Hum. Mol. Genet.* 9: 2523-2530.
6. Coenen, M.J., Smeitink, J.A., Smeets, R., Trijbels, F.J. and van den Heuvel, L.P. 2005. Mutation detection in four candidate genes (OXA1L, MRS2L, YME1L and MIPEP) for combined deficiencies in the oxidative phosphorylation system. *J. Inher. Metab. Dis.* 28: 1091-1097.
7. Online Mendelian Inheritance in Man, OMIM™. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 229300. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: MIPEP (human) mapping to 13q12.12; Mipep (mouse) mapping to 14 D1.

SOURCE

MIPEP (S-18) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of MIPEP of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-84203 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MIPEP (S-18) is recommended for detection of MIPEP of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MIPEP (S-18) is also recommended for detection of MIPEP in additional species, including canine and porcine.

Suitable for use as control antibody for MIPEP siRNA (h): sc-106225, MIPEP siRNA (m): sc-149438, MIPEP shRNA Plasmid (h): sc-106225-SH, MIPEP shRNA Plasmid (m): sc-149438-SH, MIPEP shRNA (h) Lentiviral Particles: sc-106225-V and MIPEP shRNA (m) Lentiviral Particles: sc-149438-V.

Molecular Weight of MIPEP: 81 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.