

# ME2 (G-17): sc-85051

## BACKGROUND

ME2, malic enzyme 2, is a mitochondrial NAD-dependent malic enzyme. Found as a homotetramer, the main function of ME2 is to catalyze the anaplerotic reaction yielding pyruvate from (s)-malate. One single-nucleotide polymorphism of the ME2 gene has been associated with a phenotype manifested as psychosis. ME2 directly interacts with the malate shuttle system and has roles in neuronal synthesis of glutamate and  $\gamma$ -amino butyric acid. This system has been shown to be altered in schizophrenia and bipolar disorder, which is also characterized by forms of epilepsy brought on by the role ME2 on ion channels in the brain.

## REFERENCES

1. Stöhlmacher, P. and Haferland, W. 1989. Human malic enzyme-2 polymorphism in the GDR. *Hum. Hered.* 39: 58-60.
2. Gardiner, M. 2005. Genetics of idiopathic generalized epilepsies. *Epilepsia* 46: 15-20.
3. Lenzen, K.P., Heils, A., Lorenz, S., Hempelmann, A. and Sander, T. 2005. Association analysis of malic enzyme 2 gene polymorphisms with idiopathic generalized epilepsy. *Epilepsia* 46: 1637-1641.
4. Turnbull, J., Lohi, H., Kearney, J.A., Rouleau, G.A., Delgado-Escueta, A.V., Meisler, M.H., Cossette, P. and Minassian, B.A. 2005. Sacred disease secrets revealed: the genetics of human epilepsy. *Hum. Mol. Genet.* 14: 2491-2500.
5. Pongratz, R.L., Kibbey, R.G., Shulman, G.I. and Cline, G.W. 2007. Cytosolic and mitochondrial malic enzyme isoforms differentially control Insulin secretion. *J. Biol. Chem.* 282: 200-207.
6. Lee, B.D., Walss-Bass, C., Thompson, P.M., Dassori, A., Montero, P.A., Medina, R., Contreras, S., Armas, R., Ramirez, M., Pereira, M., Salazar, R., Leach, R.J., Quezada, P., Raventos, H. and Escamilla, M.A. 2007. Malic enzyme 2 and susceptibility to psychosis and mania. *Psychiatry Res.* 150: 1-11.
7. Tronconi, M.A., Fahnenstich, H., Gerrard Weehler, M.C., Andreo, C.S., Flügge, U.I., Drincovich, M.F. and Maurino, V.G. 2008. *Arabidopsis* NAD-malic enzyme functions as a homodimer and heterodimer and has a major impact on nocturnal metabolism. *Plant Physiol.* 146: 1540-1552.

## CHROMOSOMAL LOCATION

Genetic locus: ME2 (human) mapping to 18q21.2; Me2 (mouse) mapping to 18 E2.

## SOURCE

ME2 (G-17) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of ME2 of human origin.

## PRODUCT

Each vial contains 100  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85051 P, (100  $\mu$ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

## APPLICATIONS

ME2 (G-17) is recommended for detection of ME2 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

ME2 (G-17) is also recommended for detection of ME2 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for ME2 siRNA (h): sc-75764, ME2 siRNA (m): sc-149343, ME2 shRNA Plasmid (h): sc-75764-SH, ME2 shRNA Plasmid (m): sc-149343-SH, ME2 shRNA (h) Lentiviral Particles: sc-75764-V and ME2 shRNA (m) Lentiviral Particles: sc-149343-V.

Molecular Weight of ME2: 65 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.


 MONOS  
Satisfaction  
Guaranteed

Try **ME2 (F-5): sc-514850**, our highly recommended monoclonal alternative to ME2 (G-17).