

Rab 12 (P-17): sc-85110

BACKGROUND

The Ras-related superfamily of guanine nucleotide binding proteins includes the R-Ras, Rap, Ral/Rec and Rho/Rab subfamilies. Increasing data suggests an important role for Rab proteins in either endocytosis or in biosynthetic protein transport. The process of transporting newly synthesized proteins from the endoplasmic reticulum to various stacks of the Golgi complex and to secretory vesicles involves the movement of carrier vesicles and requires Rab protein function. Rab proteins are also an integral part of endocytic pathways. Rab 12 is a 244 amino acid protein that is anchored to the membrane of the golgi apparatus and belongs to the Rab family of GTPase proteins. Like other Rab proteins, Rab 12 is thought to play a role in protein transport and may participate in vesicular trafficking events.

REFERENCES

1. Elferink, L.A., Anzai, K. and Scheller, R.H. 1992. rab15, a novel low molecular weight GTP-binding protein specifically expressed in rat brain. *J. Biol. Chem.* 267: 5768-5775.
2. Olkkonen, V.M., Dupree, P., Killisch, I., Lütcke, A., Zerial, M. and Simons, K. 1993. Molecular cloning and subcellular localization of three GTP-binding proteins of the rab subfamily. *J. Cell Sci.* 106: 1249-1261.
3. Iida, H., Wang, L., Nishii, K., Ookuma, A. and Shibata, Y. 1996. Identification of rab12 as a secretory granule-associated small GTP-binding protein in atrial myocytes. *Circ. Res.* 78: 343-347.
4. Ishido, M. and Masuo, Y. 2004. Transcriptome of pituitary adenylate cyclase-activating polypeptide-differentiated PC-12 cells. *Regul. Pept.* 123: 15-21.
5. Iida, H., Noda, M., Kaneko, T., Doiguchi, M. and Mori, T. 2005. Identification of rab12 as a vesicle-associated small GTPase highly expressed in Sertoli cells of rat testis. *Mol. Reprod. Dev.* 71: 178-185.
6. Chakrabarty, K. and Heumann, R. 2008. Prospective of Ras signaling in stem cells. *Biol. Chem.* 389: 791-798.
7. Carriere, A., Ray, H., Blenis, J. and Roux, P.P. 2008. The RSK factors of activating the Ras/MAPK signaling cascade. *Front. Biosci.* 13: 4258-4275.

CHROMOSOMAL LOCATION

Genetic locus: RAB12 (human) mapping to 18p11.22.

SOURCE

Rab 12 (P-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of Rab 12 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85110 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Rab 12 (P-17) is recommended for detection of Rab 12 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other Rab family members .

Suitable for use as control antibody for Rab 12 siRNA (h): sc-76311, Rab 12 shRNA Plasmid (h): sc-76311-SH and Rab 12 shRNA (h) Lentiviral Particles: sc-76311-V.

Molecular Weight of Rab 12: 27 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.