

ST8Sia III (S-17): sc-85169

BACKGROUND

ST8Sia III (ST8 α -N-acetyl-neuraminide α -2,8-sialyltransferase 3), also known as SIAT8C, is a 380 amino acid single-pass type II membrane protein that localizes to the membrane of the Golgi apparatus. Expressed in fetal liver and fetal and adult brain, the expression of ST8Sia III is tissue-specific and developmentally regulated. ST8Sia III plays an important role in protein modification and glycosylation, and functions to catalyze the transfer of sialic acid through α -2,8-linkage to intact fetuin glycoprotein. ST8Sia III can form polysialic acid (PSA) *in vitro* directly on α -2,3-, α -2,6- or α -2,8-linked sialic acid. PSA is an important regulator of neuronal plasticity and is present in embryonic brain tissue, where it interacts with NCAM (neural cell adhesion molecule) and plays a crucial role in fetal brain development. The gene encoding ST8Sia III maps to human chromosome 18q21.31.

REFERENCES

- Zeng, G., et al. 1997. Cloning of the cDNA coding for rat brain CMP-NeuAc:GD3 α 2-8 sialyltransferase. *Gene* 187: 131-134.
- Lee, Y.C., et al. 1998. Cloning and expression of cDNA for a human Sia α 2,3Gal β 1, 4GlcNA: α 2,8-sialyltransferase (hST8Sia III). *Arch. Biochem. Biophys.* 360: 41-46.
- Angata, K., et al. 2000. Differential biosynthesis of polysialic acid on neural cell adhesion molecule (NCAM) and oligosaccharide acceptors by three distinct α 2,8-sialyltransferases, ST8Sia IV (PST), ST8Sia II (STX), and ST8Sia III. *J. Biol. Chem.* 275: 18594-18601.
- Angata, K., et al. 2003. Polysialyltransferases: major players in polysialic acid synthesis on the neural cell adhesion molecule. *Biochimie* 85: 195-206.
- Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 609478. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Kim, S.J., et al. 2006. Molecular mechanisms involved in transcriptional activation of the human Sia- α 2,3-Gal- β 1,4-GlcNAc-R: α 2,8-sialyltransferase (hST8Sia III) gene induced by KCl in human glioblastoma cells. *Biochem. Biophys. Res. Commun.* 344: 1057-1064.
- Bork, K., et al. 2007. Experimental approaches to interfere with the polysialylation of the neural cell adhesion molecule *in vitro* and *in vivo*. *J. Neurochem.* 103: 65-71.

CHROMOSOMAL LOCATION

Genetic locus: ST8SIA3 (human) mapping to 18q21.31; St8sia3 (mouse) mapping to 18 E1.

SOURCE

ST8Sia III (S-17) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of ST8Sia III of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85169 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

ST8Sia III (S-17) is recommended for detection of ST8Sia III of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with ST8Sia II, ST8Sia IV, and ST8Sia V.

ST8Sia III (S-17) is also recommended for detection of ST8Sia III in additional species, including canine and avian.

Suitable for use as control antibody for ST8Sia III siRNA (h): sc-76580, ST8Sia III siRNA (m): sc-153867, ST8Sia III shRNA Plasmid (h): sc-76580-SH, ST8Sia III shRNA Plasmid (m): sc-153867-SH, ST8Sia III shRNA (h) Lentiviral Particles: sc-76580-V and ST8Sia III shRNA (m) Lentiviral Particles: sc-153867-V.

Molecular Weight of ST8Sia III: 44 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.