



cystatin 11 (T-19): sc-85517

BACKGROUND

The cystatins comprise a superfamily of proteins that contain multiple cystatin-like sequences and may function as active cysteine protease inhibitors. Cystatin proteins are divided into three distinct subfamilies, designated type 1 cystatins (stefins), type 2 cystatins and kininogens. The type 2 cystatins are expressed in human fluids and secretions and are encoded by genes that localize to a cystatin locus on chromosome 20. Cystatin 11, also known as SC13 or CST8L, is a 138 amino acid protein that localizes to both the nucleus and the cytoplasm and belongs to the type 2 subfamily of cystatins. Detected in epididymal and epithelial tissue, cystatin 11 functions as an antibacterial protein that is thought to play a role in sperm maturation and fertilization. Multiple isoforms of cystatin 11 exist due to alternative splicing events.

REFERENCES

- Saitoh, E., Isemura, S., Sanada, K., Kim, H.S., Smithies, O. and Maeda, N. 1988. Cystatin superfamily. Evidence that family II cystatin genes are evolutionarily related to family III cystatin genes. *Biol. Chem. Hoppe-Seyler.* 369 Suppl: 191-197.
- Dickinson, D.P., Thiesse, M., Dempsey, L.D. and Millar, S.J. 1993. Genomic cloning, physical mapping, and expression of human type 2 cystatin genes. *Crit. Rev. Oral Biol. Med.* 4: 573-580.
- Thiesse, M., Millar, S.J. and Dickinson, D.P. 1994. The human type 2 cystatin gene family consists of eight to nine members, with at least seven genes clustered at a single locus on human chromosome 20. *DNA Cell Biol.* 13: 97-116.
- Dickinson, D.P., Zhao, Y., Thiesse, M. and Siciliano, M.J. 1994. Direct mapping of seven genes encoding human type 2 cystatins to a single site located at 20p11.2. *Genomics* 24: 172-175.
- Brown, W.M. and Dziejglewska, K.M. 1997. Friends and relations of the cystatin superfamily—new members and their evolution. *Protein Sci.* 6: 5-12.
- Hamil, K.G., Liu, Q., Sivashanmugam, P., Yenugu, S., Soundararajan, R., Grossman, G., Richardson, R.T., Zhang, Y.L., O'Rand, M.G., Petrusz, P., French, F.S. and Hall, S.H. 2002. Cystatin 11: a new member of the cystatin type 2 family. *Endocrinology* 143: 2787-2796.
- Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609731: World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Cornwall, G.A. and Hsia, N. 2003. A new subgroup of the family 2 cystatins. *Mol. Cell. Endocrinol.* 200: 1-8.
- Frygeliuss, J., Oscarson, M., Nordqvist, K., Wedell, A. and Töhönen, V. 2007. The reproductive tissue specific cystatin subgroup of genes: expression during gonadal development in wildtype and testatin knockout animals. *Sex Dev.* 1: 363-372.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: CST11 (human) mapping to 20p11.21.

SOURCE

cystatin 11 (T-19) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of cystatin 11 of human origin.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85517 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

cystatin 11 (T-19) is recommended for detection of cystatin 11 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with isoform cystatin-2.

Suitable for use as control antibody for cystatin 11 siRNA (h): sc-77080, cystatin 11 shRNA Plasmid (h): sc-77080-SH and cystatin 11 shRNA (h) Lentiviral Particles: sc-77080-V.

Molecular Weight of cystatin 11: 17 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.