

LLC1 (P-16): sc-85607

BACKGROUND

LLC1 (low in lung cancer 1), also known as C20orf85, is a 137 amino acid protein that may play a role in lung tumorigenesis. The gene encoding LLC1 is located on human chromosome 20, which is comprised of approximately 2% of the human genome, chromosome 20 contains nearly 63 million bases that encode over 600 genes, some of which are associated with Creutzfeldt-Jakob disease, amyotrophic lateral sclerosis, spinal muscular atrophy, ring chromosome 20 epilepsy syndrome and Alagille syndrome. Additionally, chromosome 20 contains a region with numerous genes which are thought important for seminal production and may be potential targets for male contraception.

REFERENCES

- Joó, J.G., Beke, A., Tóth-Pál, E., Hargitai, B., Szigeti, Z., Papp, C. and Papp, Z. 2006. Trisomy 20 mosaicism and nonmosaic trisomy 20: a report of 2 cases. *J. Reprod. Med.* 51: 209-212.
- Ville, D., Kaminska, A., Bahi-Buisson, N., Biraben, A., Plouin, P., Telvi, L., Dulac, O. and Chiron, C. 2006. Early pattern of epilepsy in the ring chromosome 20 syndrome. *Epilepsia* 47: 543-549.
- Elghezal, H., Hannachi, H., Mougou, S., Kammoun, H., Triki, C. and Saad, A. 2007. Ring chromosome 20 syndrome without deletions of the subtelomeric and CHRNA4-KCNQ2 genes loci. *Eur. J. Med. Genet.* 50: 441-445.
- Kazantsev, A.G. 2007. Cellular pathways leading to neuronal dysfunction and degeneration. *Drug News Perspect.* 20: 501-509.
- Lundwall, A. 2007. A locus on chromosome 20 encompassing genes that are highly expressed in the epididymis. *Asian J. Androl.* 9: 540-544.
- Hong, K.M., Yang, S.H., Chowdhuri, S.R., Player, A., Hames, M., Fukuoka, J., Meerzaman, D., Dracheva, T., Sun, Z., Yang, P. and Jen, J. 2007. Inactivation of LLC1 gene in nonsmall cell lung cancer. *Int. J. Cancer.* 120: 2353-2358.
- Savai, R., Langheinrich, A.C., Schermuly, R.T., Pullamsetti, S.S., Dumitrascu, R., Traupe, H., Rau, W.S., Seeger, W., Grimminger, F. and Banat, G.A. 2009. Evaluation of angiogenesis using micro-computed tomography in a xenograft mouse model of lung cancer. *Neoplasia* 11: 48-56.
- Lee, D.Y., Lee, S.W., Kim, S.K., Lee, M., Chang, H.W., Moon, H.T., Byun, Y. and Kim, S.Y. 2009. Antiangiogenic activity of orally absorbable heparin derivative in different types of cancer cells. *Pharm. Res.* 26: 2667-2676.

CHROMOSOMAL LOCATION

Genetic locus: C20orf85 (human) mapping to 20q13.32.

SOURCE

LLC1 (P-16) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of LLC1 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 100 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85607 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

LLC1 (P-16) is recommended for detection of LLC1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

LLC1 (P-16) is also recommended for detection of LLC1 in additional species, including bovine and porcine.

Suitable for use as control antibody for LLC1 siRNA (h): sc-75428, LLC1 shRNA Plasmid (h): sc-75428-SH and LLC1 shRNA (h) Lentiviral Particles: sc-75428-V.

Molecular Weight of LLC1: 16 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.