

pki γ (E-16): sc-85825

BACKGROUND

The second messenger cyclic AMP (cAMP) mediates a diverse array of cellular responses such as proliferation, ion transport, regulation of metabolism and gene transcription by activating the cAMP-dependent protein kinase (cAPK or PKA). Protein kinase inhibitors are potent inhibitors of the catalytic subunit of PKAs. Pki γ (cAMP-dependent protein kinase inhibitor γ) is a 76 amino acid protein that inhibits catalytic subunit-dependent transcription and efficiently terminates nuclear PKA activity. Knockdown of the mRNA encoding pki γ substantially extends the anti-apoptotic effect of parathyroid hormone and β -adrenergic agonists. This evidence suggests that inhibition of pki γ activity may be a useful co-therapy in the treatment of osteoporosis. Pki γ is highly expressed in skeletal muscle, testis and heart.

REFERENCES

- Collins, S.P. and Uhler, M.D. 1997. Characterization of pki γ , a novel isoform of the protein kinase inhibitor of cAMP-dependent protein kinase. *J. Biol. Chem.* 272: 18169-18178.
- Zheng, L., et al. 2000. Cloning and mapping of human pki β and pki γ , and comparison of tissue expression patterns of three members of the protein kinase inhibitor family, including pki α . *Biochem. J.* 349: 403-407.
- Taylor, S.S., et al. 2004. PKA: a portrait of protein kinase dynamics. *Biochim. Biophys. Acta* 1697: 259-269.
- Online Mendelian Inheritance in Man, OMIM™. 2004. Johns Hopkins University, Baltimore, MD. MIM Number: 604932. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Taylor, S.S., et al. 2005. Dynamics of signaling by PKA. *Biochim. Biophys. Acta* 1754: 25-37.
- Chen, X., et al. 2005. Endogenous protein kinase inhibitor γ terminates immediate-early gene expression induced by cAMP-dependent protein kinase (PKA) signaling: termination depends on PKA inactivation rather than PKA export from the nucleus. *J. Biol. Chem.* 280: 2700-2707.

CHROMOSOMAL LOCATION

Genetic locus: PKIG (human) mapping to 20q13.12; Pkig (mouse) mapping to 2 H3.

SOURCE

pki γ (E-16) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of pki γ of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-85825 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

pki γ (E-16) is recommended for detection of pki γ of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

pki γ (E-16) is also recommended for detection of pki γ in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for pki γ siRNA (h): sc-106417, pki γ siRNA (m): sc-155936, pki γ shRNA Plasmid (h): sc-106417-SH, pki γ shRNA Plasmid (m): sc-155936-SH, pki γ shRNA (h) Lentiviral Particles: sc-106417-V and pki γ shRNA (m) Lentiviral Particles: sc-155936-V.

Molecular Weight of pki γ : 8 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfation
Guaranteed

Try **pki γ (D-8): sc-515108**, our highly recommended monoclonal alternative to pki γ (E-16).