

GGTL3 (N-13): sc-86143

BACKGROUND

Glutathione is a tripeptide electron donor that functions as an antioxidant, protecting cells from toxins by reducing disulfide bonds during oxidative stress. The metabolism of glutathione requires a variety of enzymes, such as GGT1, GGT2, GGT5, GGTL2 and GGTL3. GGTL3, also known as GGT7 (γ -glutamyltransferase 7), GGT4 or GGTL5, is a 662 amino acid single-pass type II membrane protein that belongs to the γ -glutamyltransferase family and is involved in the pathway of Glutathione metabolism. Expressed at low levels in multiple tissues, including lung, liver, heart, brain, testis and spleen, GGTL3 functions as a heterodimer of light and heavy chains that work together to cleave Glutathione-conjugated peptides, thereby releasing an unconjugated protein. GGTL3 exists as multiple alternatively spliced isoforms and, via its catalytic activity, may be involved in the pathogenesis of lung cancer.

REFERENCES

1. Courtay, C., Heisterkamp, N., Siest, G. and Groffen, J. 1994. Expression of multiple γ -glutamyltransferase genes in man. *Biochem. J.* 297: 503-508.
2. He, X., Di, Y., Li, J., Xie, Y., Tang, Y., Zhang, F., Wei, L., Zhang, Y., Qin, W., Huo, K., Li, Y., Wan, D. and Gu, J. 2002. Molecular cloning and characterization of CT120, a novel membrane-associated gene involved in amino acid transport and glutathione metabolism. *Biochem. Biophys. Res. Commun.* 297: 528-536.
3. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 612342. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Llamazares, M., Cal, S., Quesada, V. and López-Otín, C. 2003. Identification and characterization of ADAMTS-20 defines a novel subfamily of metalloproteinases-disintegrins with multiple Thrombospondin 1 repeats and a unique GON domain. *J. Biol. Chem.* 278: 13382-13389.
5. Wu, G., Fang, Y.Z., Yang, S., Lupton, J.R. and Turner, N.D. 2004. Glutathione metabolism and its implications for health. *J. Nutr.* 134: 489-492.
6. Ristoff, E. and Larsson, A. 2007. Inborn errors in the metabolism of Glutathione. *Orphanet. J. Rare Dis.* 2: 16.
7. Heisterkamp, N., Groffen, J., Warburton, D. and Sneddon, T.P. 2008. The human γ -glutamyltransferase gene family. *Hum. Genet.* 123: 321-332.

CHROMOSOMAL LOCATION

Genetic locus: GGT7 (human) mapping to 20q11.22; Ggt7 (mouse) mapping to 2 H1.

SOURCE

GGTL3 (N-13) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of GGTL3 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-86143 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

GGTL3 (N-13) is recommended for detection of GGTL3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with GGTL3-4.

GGTL3 (N-13) is also recommended for detection of GGTL3 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for GGTL3 siRNA (h): sc-75129, GGTL3 siRNA (m): sc-145393, GGTL3 shRNA Plasmid (h): sc-75129-SH, GGTL3 shRNA Plasmid (m): sc-145393-SH, GGTL3 shRNA (h) Lentiviral Particles: sc-75129-V and GGTL3 shRNA (m) Lentiviral Particles: sc-145393-V.

Molecular Weight of GGTL3: 70 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.