

CAT-4 (C-15): sc-86351

BACKGROUND

As a member of the APC family of transporters, CAT-4 (cationic amino acid transporter 4), also known as solute carrier family 7 member 4, is a 635 amino acid multi-pass membrane protein that is involved in the transport of the cationic amino acids, Arginine, Lysine and Ornithine. This uptake of cationic amino acids is designated "system Y⁺", which is pH insensitive, stereoselective, Na⁺ independent and inhibitable by neural amino acids in the presence of Na⁺. CAT-4 displays high sequence similarity with CAT-1 and CAT-2 and is expressed in testis, placenta and brain. A microdeletion of the chromosomal band near the location of the gene encoding CAT-4 causes velocardiofacial syndrome (also known as DiGeorge syndrome), a disease that is characterized by several clinical findings, including conotruncal cardiac defects, prominent tubular nose and a hypernasal voice. This suggests that CAT-4 may play a role in determining the velocardiofacial phenotype.

REFERENCES

1. Sperandio, M.P., Borsani, G., Incerti, B., Zollo, M., Rossi, E., Zuffardi, O., Castaldo, P., Tagliatela, M., Andria, G. and Sebastio, G. 1998. The gene encoding a cationic amino acid transporter (SLC7A4) maps to the region deleted in the velocardiofacial syndrome. *Genomics* 49: 230-236.
2. Hammermann, R., Brunn, G. and Racke, K. 2001. Analysis of the genomic organization of the human cationic amino acid transporters CAT-1, CAT-2 and CAT-4. *Amino Acids*. 21: 211-219.
3. Wolf, S., Janzen, A., Vekony, N., Martine, U., Strand, D. and Closs, E.I. 2002. Expression of solute carrier 7A4 (SLC7A4) in the plasma membrane is not sufficient to mediate amino acid transport activity. *Biochem. J.* 364: 767-775.
4. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 192430: World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Rotoli, B.M., Bussolati, O., Sala, R., Gazzola, G.C. and Dall'Asta, V. 2005. The transport of cationic amino acids in human airway cells: expression of system Y⁺L activity and transepithelial delivery of NOS inhibitors. *FASEB J.* 19: 810-812.
6. Closs, E.I., Boissel, J.P., Habermeier, A. and Rotmann, A. 2006. Structure and function of cationic amino acid transporters (CATs). *J. Membr. Biol.* 213: 67-77.
7. Emanuel, B.S. 2008. Molecular mechanisms and diagnosis of chromosome 22q11.2 rearrangements. *Dev. Disabil. Res. Rev.* 14: 11-18.

CHROMOSOMAL LOCATION

Genetic locus: SLC7A4 (human) mapping to 22q11.21.

SOURCE

CAT-4 (C-15) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of CAT-4 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-86351 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CAT-4 (C-15) is recommended for detection of CAT-4 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with family members CAT-1 or CAT-2.

Suitable for use as control antibody for CAT-4 siRNA (h): sc-72804, CAT-4 shRNA Plasmid (h): sc-72804-SH and CAT-4 shRNA (h) Lentiviral Particles: sc-72804-V.

Molecular Weight of CAT-4: 68 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker[™] compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker[™] Molecular Weight Standards: sc-2035, TBS Blotting A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz[™] Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfaction
Guaranteed

Try **CAT-4 (B-2): sc-376557**, our highly recommended monoclonal alternative to CAT-4 (C-15).