

DDX17 (N-13): sc-86408

-BACKGROUND

DEAD box proteins, characterized by the conserved motif ASP-Glu-Ala-ASP, are putative RNA helicases implicated in several cellular processes involving modifications of RNA secondary structure. Specifically, DEAD box proteins are involved in translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, members of this family may be involved in embryogenesis, spermatogenesis and cellular growth and division. DDX17, also designated p72, is highly homologous to DDX5 (p68). DDX17 and DDX5 have been implicated in growth regulation by acting as transcriptional co-regulators for several transcription factors, including ER α , p53, MyoD and RUNX2. Impairment of DDX17 may affect early brain development that can lead to Down Syndrome. Alternatively, up-regulation of DDX17 and DDX5 directly contributes to colon cancer implicating DDX17 as a therapeutic target to combat against colon cancer.

CHROMOSOMAL LOCATION

Genetic locus: DDX17 (human) mapping to 22q13.1; Ddx17 (mouse) mapping to 15 E1.

SOURCE

DDX17 (N-13) is an affinity purified rabbit polyclonal antibody raised against a peptide mapping within an internal region of DDX17 of human origin.

PRODUCT

Each vial contains 100 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-86408 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

DDX17 (N-13) is recommended for detection of DDX17 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other DDX family members.

DDX17 (N-13) is also recommended for detection of DDX17 in additional species, including canine.

Suitable for use as control antibody for DDX17 siRNA (h): sc-77106, DDX17 siRNA (m): sc-142922, DDX17 shRNA Plasmid (h): sc-77106-SH, DDX17 shRNA Plasmid (m): sc-142922-SH, DDX17 shRNA (h) Lentiviral Particles: sc-77106-V and DDX17 shRNA (m) Lentiviral Particles: sc-142922-V.


Molecular Weight of DDX17 isoforms: 72/82 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204, NIH/3T3 whole cell lysate: sc-2210 or HEK293 whole cell lysate: sc-45136.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



DDX17 (N-13): sc-86408. Western blot analysis of DDX17 expression in Jurkat whole cell lysate.



DDX17 (N-13): sc-86408. Immunofluorescence staining of formalin-fixed HepG2 cells showing nuclear localization.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **DDX17 (H-7): sc-398168** or **DDX17 (C-9): sc-271112**, our highly recommended monoclonal alternatives to DDX17 (N-13).