



DJBP (C-17): sc-86418

BACKGROUND

DJ-1 was first identified as a novel oncogene that transforms mouse NIH/3T3 cells in cooperation with activated Ras. Also, DJ-1 positively regulates androgen receptor (AR) by forming a complex with PIASx, a negative regulator of AR. Mutations in the DJ-1 gene have been implicated in Parkinson's disease, and loss of DJ-1 function leads to neurodegeneration. DJBP (DJ-1-binding protein), also known as EF-hand calcium-binding domain-containing protein 6, is a 1,501 amino acid protein that binds the C-terminal region of DJ-1, as well as the DNA binding region of AR. With specific expression in testis, DJBP forms a ternary complex with DJ-1 and AR and colocalizes to the nucleus. Unlike PIASx, DJBP1 acts as an antagonist of DJ-1 by negatively regulating testosterone-dependent AR transactivation activity via recruitment of a histone deacetylase complex. There are six isoforms of DJBP which exist as a result of alternative splicing events.

REFERENCES

- Hirosawa, M., Nagase, T., Murahashi, Y., Kikuno, R. and Ohara, O. 2001. Identification of novel transcribed sequences on human chromosome 22 by expressed sequence tag mapping. *DNA Res.* 8: 1-9.
- Takahashi, K., Taira, T., Niki, T., Seino, C., Iguchi-Arigo, S.M. and Ariga, H. 2001. DJ-1 positively regulates the androgen receptor by impairing the binding of PIASx alpha to the receptor. *J. Biol. Chem.* 276: 37556-37563.
- Huai, Q., Sun, Y., Wang, H., Chin, L.S., Li, L., Robinson, H. and Ke, H. 2003. Crystal structure of DJ-1/RS and implication on familial Parkinson's disease. *FEBS Lett.* 549: 171-175.
- Niki, T., Takahashi-Niki, K., Taira, T., Iguchi-Arigo, S.M. and Ariga, H. 2003. DJBP: a novel DJ-1-binding protein, negatively regulates the androgen receptor by recruiting histone deacetylase complex, and DJ-1 antagonizes this inhibition by abrogation of this complex. *Mol. Cancer Res.* 1: 247-261.
- Taira, T., Iguchi-Arigo, S.M. and Ariga, H. 2004. Co-localization with DJ-1 is essential for the androgen receptor to exert its transcription activity that has been impaired by androgen antagonists. *Biol. Pharm. Bull.* 27: 574-577.
- Hong, C.Y., Gong, E.Y., Kim, K., Suh, J.H., Ko, H.M., Lee, H.J., Choi, H.S. and Lee, K. 2005. Modulation of the expression and transactivation of androgen receptor by the basic helix-loop-helix transcription factor Pod-1 through recruitment of histone deacetylase 1. *Mol. Endocrinol.* 19: 2245-2257.
- Tillman, J.E., Yuan, J., Gu, G., Fazli, L., Ghosh, R., Flynt, A.S., Gleave, M., Rennie, P.S. and Kasper, S. 2007. DJ-1 binds androgen receptor directly and mediates its activity in hormonally treated prostate cancer cells. *Cancer Res.* 67: 4630-4637.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: EFCAB6 (human) mapping to 22q13.2; Efcab6 (mouse) mapping to 15 E2.

SOURCE

DJBP (C-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of DJBP of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-86418 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

DJBP (C-17) is recommended for detection of DJBP of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with DJBP-5.

Suitable for use as control antibody for DJBP siRNA (h): sc-77155, DJBP siRNA (m): sc-143052, DJBP shRNA Plasmid (h): sc-77155-SH, DJBP shRNA Plasmid (m): sc-143052-SH, DJBP shRNA (h) Lentiviral Particles: sc-77155-V and DJBP shRNA (m) Lentiviral Particles: sc-143052-V.

Molecular Weight of DJBP: 173 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.