

MEI1 (D-14): sc-86731

BACKGROUND

The predominant cause of spermatogenic arrest of meiosis is the failure of homologous chromosomes to accurately synapse. MEI1 (meiosis inhibitor protein 1), also designated meiosis defective protein 1, is a 1,274 amino acid protein that is likely required for the formation of genetically programmed double-strand breaks, the first step in the initiation of meiosis. With predominant expression in testes, it is likely that defects of the gene encoding MEI1 results in male infertility. Interestingly, studies show that genetic variation in the MEI gene possibly predisposes European Americans but not Israeli men to infertility by meiotic arrest. Human MEI1 shares 79% sequence similarity with its mouse homolog. There are seven isoforms of MEI1 that are produced as a result of alternative splicing events.

REFERENCES

- Libby, B.J., De La Fuente, R., O'Brien, M.J., Wigglesworth, K., Cobb, J., Inselman, A., Eaker, S., Handel, M.A., Eppig, J.J. and Schimenti, J.C. 2002. The mouse meiotic mutation Mei1 disrupts chromosome synapsis with sexually dimorphic consequences for meiotic progression. *Dev. Biol.* 242: 174-187.
- Grelon, M., Gendrot, G., Vezon, D., Pelletier, G., Mathilde, G., Ghislaine, G., Daniel, V. and Georges, P. 2003. The *Arabidopsis* MEI1 gene encodes a protein with five BRCT domains that is involved in meiosis-specific DNA repair events independent of SPO11-induced DSBs. *Plant J.* 35: 465-475.
- Libby, B.J., Reinholdt, L.G. and Schimenti, J.C. 2003. Positional cloning and characterization of Mei1, a vertebrate-specific gene required for normal meiotic chromosome synapsis in mice. *Proc. Natl. Acad. Sci. USA* 100: 15706-15711.
- Reinholdt, L.G. and Schimenti, J.C. 2005. Mei1 is epistatic to Dmc1 during mouse meiosis. *Chromosoma* 114: 127-134.
- Liebe, B., Petukhova, G., Barchi, M., Bellani, M., Braselmann, H., Nakano, T., Pandita, T.K., Jasin, M., Fornace, A., Meistrich, M.L., Baarends, W.M., Schimenti, J., de Lange, T., Keeney, S., Camerini-Otero, R.D. and Scherthan, H. 2006. Mutations that affect meiosis in male mice influence the dynamics of the mid-preleptotene and bouquet stages. *Exp. Cell Res.* 312: 3768-3781.
- Sato, H., Miyamoto, T., Yogev, L., Namiki, M., Koh, E., Hayashi, H., Sasaki, Y., Ishikawa, M., Lamb, D.J., Matsumoto, N., Birk, O.S., Niikawa, N. and Sengoku, K. 2006. Polymorphic alleles of the human MEI1 gene are associated with human azoospermia by meiotic arrest. *J. Hum. Genet.* 51: 533-540.

CHROMOSOMAL LOCATION

Genetic locus: MEI1 (human) mapping to 22q13.2; Mei1 (mouse) mapping to 15 E1.

SOURCE

MEI1 (D-14) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of MEI1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-86731 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

MEI1 (D-14) is recommended for detection of MEI1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

MEI1 (D-14) is also recommended for detection of MEI1 in additional species, including porcine.

Suitable for use as control antibody for MEI1 siRNA (h): sc-75770, MEI1 siRNA (m): sc-149362, MEI1 shRNA Plasmid (h): sc-75770-SH, MEI1 shRNA Plasmid (m): sc-149362-SH, MEI1 shRNA (h) Lentiviral Particles: sc-75770-V and MEI1 shRNA (m) Lentiviral Particles: sc-149362-V.

Molecular Weight of MEI1: 141 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfation
Guaranteed

Try **MEI1 (F-4): sc-515359**, our highly recommended monoclonal alternative to MEI1 (D-14).