

FBXO16 (C-17): sc-87120

BACKGROUND

Belonging to the F-box family of proteins, FBXO16 (F-box only protein 16), also known as FBX16, is a 292 amino acid protein that contains one C-terminal F-box domain. F-box proteins are critical components of the SCF (Skp1-CUL-1-F-box protein)-type E3 ubiquitin ligase complex and are involved in substrate recognition and recruitment for ubiquitination. They are members of a larger family of proteins that are involved in the regulation of a wide variety of cellular mechanisms, including the cell cycle, the immune response, signaling cascades and developmental processes. They function by targeting proteins, such as cyclins, cyclin-dependent kinase inhibitors, $\text{I}\kappa\text{B-}\alpha$ and β -catenin, for degradation by the proteasome after ubiquitination. Via its F-box domain, FBXO16 can directly interact with Skp1 p19 and CUL-1. FBXO16 is expressed in heart, spleen and colon.

REFERENCES

- Cenciarelli, C., Chiaur, D.S., Guardavaccaro, D., Parks, W., Vidal, M. and Pagano, M. 1999. Identification of a family of human F-box proteins. *Curr. Biol.* 9: 1177-1179.
- Winston, J.T., Koepf, D.M., Zhu, C., Elledge, S.J. and Harper, J.W. 1999. A family of mammalian F-box proteins. *Curr. Biol.* 9: 1180-1182.
- Winston, J.T., Strack, P., Beer-Romero, P., Chu, C.Y., Elledge, S.J. and Harper, J.W. 1999. The SCFb-TRCP-ubiquitin ligase complex associates specifically with phosphorylated destruction motifs in $\text{I}\kappa\text{B-}\alpha$ and β -catenin and stimulates $\text{I}\kappa\text{B-}\alpha$ ubiquitination *in vitro*. *Genes Dev.* 13: 270-283.
- Craig, K.L. and Tyers, M. 1999. The F-box: a new motif for ubiquitin dependent proteolysis in cell cycle regulation and signal transduction. *Prog. Biophys. Mol. Biol.* 72: 299-328.
- Ilyin, G.P., Rialland, M., Pigeon, C. and Guguen-Guillouzo, C. 2000. cDNA cloning and expression analysis of new members of the mammalian F-box protein family. *Genomics* 67: 40-47.
- Schulman, B.A., Carrano, A.C., Jeffrey, P.D., Bowen, Z., Kinnucan, E.R., Finnin, M.S., Elledge, S.J., Harper, J.W., Pagano, M. and Pavletich, N.P. 2000. Insights into SCF ubiquitin ligases from the structure of the Skp1-Skp2 complex. *Nature* 408: 381-386.
- Ilyin, G.P., Serandour, A.L., Pigeon, C., Rialland, M., Glaise, D. and Guguen-Guillouzo, C. 2002. A new subfamily of structurally related human F-box proteins. *Gene* 296: 11-20.
- Jin, J., Cardozo, T., Lovering, R.C., Elledge, S.J., Pagano, M. and Harper, J.W. 2004. Systematic analysis and nomenclature of mammalian F-box proteins. *Genes Dev.* 18: 2573-2580.
- Okita, K., Ichisaka, T. and Yamanaka, S. 2007. Generation of germline-competent induced pluripotent stem cells. *Nature* 448: 313-317.

CHROMOSOMAL LOCATION

Genetic locus: FBXO16 (human) mapping to 8p21.1; Fbxo16 (mouse) mapping to 14 D1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

SOURCE

FBXO16 (C-17) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of FBXO16 of human origin.

PRODUCT

Each vial contains 200 μg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-87120 P, (100 μg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

FBXO16 (C-17) is recommended for detection of FBXO16 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μg per 100-500 μg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000); non cross-reactive with other FBXO family members.


FBXO16 (C-17) is also recommended for detection of FBXO16 in additional species, including canine and porcine.

Suitable for use as control antibody for FBXO16 siRNA (h): sc-77769, FBXO16 siRNA (m): sc-145106, FBXO16 shRNA Plasmid (h): sc-77769-SH, FBXO16 shRNA Plasmid (m): sc-145106-SH, FBXO16 shRNA (h) Lentiviral Particles: sc-77769-V and FBXO16 shRNA (m) Lentiviral Particles: sc-145106-V.

Molecular Weight of FBXO16: 35 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200, NIH/3T3 whole cell lysate: sc-2210 or SW480 cell lysate: sc-2219.

DATA



FBXO16 (C-17): sc-87120. Western blot analysis of FBXO16 expression in SW480 (A), HeLa (B) and NIH/3T3 (C) whole cell lysates.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.