

CD98 (H-300): sc-9160

BACKGROUND

CD98 (4F2, CD98, MDU1, 4F2HC, 4T2HC, NACAE) is a disulfide-linked heterodimer composed of a glycosylated heavy chain and a non-glycosylated light chain. CD98 is a scaffolding protein that interacts with basolaterally expressed amino acid transporters and β 1 integrins and can alter amino acid transport and cell adhesion, migration and branching morphogenesis. The heavy chain is a type II integral membrane protein. CD98 is expressed on T cells and is upregulated upon T cell activation. CD98 is also present on monocytes and at lower levels on granulocytes, platelets and lymphocytes. Evidence suggests that CD98 may play a role in the regulation of T cell activation and proliferation. Alternate transcriptional splice variants, encoding different isoforms, exist for the human CD98 gene.

CHROMOSOMAL LOCATION

Genetic locus: SLC3A2 (human) mapping to 11q12.3; Slc3a2 (mouse) mapping to 19 A.

SOURCE

CD98 (H-300) is a rabbit polyclonal antibody raised against amino acids 230-529 of CD98 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

CD98 (H-300) is recommended for detection of CD98 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CD98 siRNA (h): sc-35033, CD98 siRNA (m): sc-35034, CD98 shRNA Plasmid (h): sc-35033-SH, CD98 shRNA Plasmid (m): sc-35034-SH, CD98 shRNA (h) Lentiviral Particles: sc-35033-V and CD98 shRNA (m) Lentiviral Particles: sc-35034-V.

Molecular Weight of CD98: 125 kDa.

Positive Controls: Hep G2 cell lysate: sc-2227, HeLa whole cell lysate: sc-2200 or U-937 cell lysate: sc-2239.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

DATA



CD98 (H-300): sc-9160. Western blot analysis of CD98 expression in U-937 (A), HeLa (B) and Hep G2 (C) whole cell lysates.

CD98 (H-300): sc-9160. Immunofluorescence staining of methanol-fixed U-937 cells showing membrane localization (A). Immunoperoxidase staining of formalin fixed, paraffin-embedded human parathyroid gland tissue showing membrane staining of glandular cells magnification. Kindly provided by The Swedish Human Protein Atlas (HPA) program (B).

SELECT PRODUCT CITATIONS

- Kudo, Y., et al. 2003. Manipulation of CD98 expression affects both trophoblast cell fusion and amino acid transport activity during syncytialization of human placental BeWo cells. *J. Physiol.* 550: 3-9.
- Lemaitre, G., et al. 2005. CD98, a novel marker of transient amplifying human keratinocytes. *Proteomics* 5: 3637-3645.
- Kierszenbaum, A.L., et al. 2006. Role of integrins, tetraspanins, and ADAM proteins during the development of apoptotic bodies by spermatogenic cells. *Mol. Reprod. Dev.* 73: 906-917.
- Dalton, P., et al. 2007. Membrane trafficking of CD98 and its ligand galectin 3 in BeWo cells—implication for placental cell fusion. *FEBS J.* 274: 2715-2727.
- Dalton, P., et al. 2007. Differential effect of cross-linking the CD98 heavy chain on fusion and amino acid transport in the human placental trophoblast (BeWo) cell line. *Biochim. Biophys. Acta* 1768: 401-410.
- Ramya, T.N., et al. 2010. *In situ trans* ligands of CD22 identified by glycan-protein photocross-linking-enabled proteomics. *Mol. Cell. Proteomics* 9: 1339-1351.
- Kumar, R., et al. 2010. Metastasis-associated protein 1 short form stimulates Wnt1 pathway in mammary epithelial and cancer cells. *Cancer Res.* 70: 6598-6608.
- Isono, T. 2011. O-GlcNAc-specific antibody CTD110.6 cross-reacts with N-GlcNAc2-modified proteins induced under glucose deprivation. *PLoS ONE* 6: e18959.
- Liu, Y., et al. 2011. An immunohistochemical analysis-based decision tree model for estimating the risk of lymphatic metastasis in pN0 squamous cell carcinomas of the lung. *Histopathology* 59: 882-891.
- Isono, T., et al. 2013. Study of global transcriptional changes of N-GlcNAc2 proteins-producing T24 bladder carcinoma cells under glucose deprivation. *PLoS ONE* 8: e60397.