

# INSM2 siRNA (h): sc-92409

## BACKGROUND

Insulinoma-associated protein 2 (INSM2), also known as zinc-finger protein IA-6, is a 566 amino acid nuclear protein. INSM2 is believed to act as a growth or tumor suppressor in certain neurons and in liver cells. INSM2 contains five zinc-finger C2H2-type domains which bind nucleic acids. Zinc-fingers are composed of 25-30 amino acid residues and have the ability to bind about five nucleotides. C2H2-type zinc-fingers have two cysteines on one extremity of the domain, and two histidines on the other extremity. These four cysteines and histidines interact with one zinc atom, resulting in the finger-like conformation of the domain. It is believed that the C2H2-type zinc-fingers require the binding of a zinc atom in order for the finger to bind DNA or RNA.

## REFERENCES

1. Pellegrino, G.R. and Berg, J.M. 1991. Identification and characterization of "zinc-finger" domains by the polymerase chain reaction. *Proc. Natl. Acad. Sci. USA* 88: 671-675.
2. Hoffman, R.C., Horvath, S.J. and Klevit, R.E. 1993. Structures of DNA-binding mutant zinc-finger domains: implications for DNA binding. *Protein Sci.* 2: 951-965.
3. Tateno, M., Fukunishi, Y., Komatsu, S., Okazaki, Y., Kawai, J., Shibata, K., Itoh, M., Muramatsu, M., Held, W.A. and Hayashizaki, Y. 2001. Identification of a novel member of the snail/Gfi-1 repressor family, mlt 1, which is methylated and silenced in liver tumors of SV40 T antigen transgenic mice. *Cancer Res.* 61: 1144-1153.
4. Shimojo, M., Lee, J.H. and Hersh, L.B. 2001. Role of zinc-finger domains of the transcription factor neuron-restrictive silencer factor/repressor element-1 silencing transcription factor in DNA binding and nuclear localization. *J. Biol. Chem.* 276: 13121-13126.
5. Cai, T., Fukushige, T., Notkins, A.L. and Krause, M. 2004. Insulinoma-associated protein IA-2, a vesicle transmembrane protein, genetically interacts with UNC-31/CAPS and affects neurosecretion in *Caenorhabditis elegans*. *J. Neurosci.* 24: 3115-3124.
6. Dreier, B., Fuller, R.P., Segal, D.J., Lund, C.V., Blancafort, P., Huber, A., Kokscho, B. and Barbas, C.F. 2005. Development of zinc-finger domains for recognition of the 5'-CNN-3' family DNA sequences and their use in the construction of artificial transcription factors. *J. Biol. Chem.* 280: 35588-35597.
7. Pedersen, N., Pedersen, M.W., Lan, M.S., Breslin, M.B. and Poulsen, H.S. 2006. The insulinoma-associated 1: a novel promoter for targeted cancer gene therapy for small-cell lung cancer. *Cancer Gene Ther.* 13: 375-384.

## CHROMOSOMAL LOCATION

Genetic locus: INSM2 (human) mapping to 14q13.2.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

INSM2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see INSM2 shRNA Plasmid (h): sc-92409-SH and INSM2 shRNA (h) Lentiviral Particles: sc-92409-V as alternate gene silencing products.

For independent verification of INSM2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-92409A, sc-92409B and sc-92409C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

INSM2 siRNA (h) is recommended for the inhibition of INSM2 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor INSM2 gene expression knockdown using RT-PCR Primer: INSM2 (h)-PR: sc-92409-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.