

MCP-4 (E-20): sc-9656

BACKGROUND

The monocyte chemotactic proteins, MCP-1, MCP-2 and MCP-3, form a sub-family of the C-C or β -chemokines which are characterized by a set of conserved adjacent cysteines. MCPs are produced by a variety of cells, including T lymphocytes, subsequent to their activation with cytokines such as IL-1, TNF α , and IFN- γ . MCP-1 levels are increased during infection and inflammation, which are both characterized by leukocyte infiltration. *In vitro* studies have shown that the MCP isoforms exhibit their chemotactic effects on different subpopulations of lymphocytes. MCP-1 is a potent basophil activator but does not effect eosinophils, whereas MCP-2 stimulates both eosinophils and basophils. MCP-3 has been shown to have the broadest range of influence, activating monocytes, dendritic cells, lymphocytes, natural killer cells, eosinophils, basophils and neutrophils. Two MCP-1 receptors that differ in their carboxy termini have been identified.

REFERENCES

1. Charo, I.F., et al. 1994. Molecular cloning and functional expression of two monocyte chemoattractant protein 1 receptors reveals alternative splicing of the carboxyl-terminal tails. *Proc. Natl. Acad. Sci. USA* 91: 2752-2756.
2. Taub, D.D., et al. 1995. Monocyte chemotactic protein-1 (MCP-1), -2, and -3 are chemotactic for human T lymphocytes. *J. Clin. Invest.* 95: 1370-1376.
3. Weber, M., et al. 1995. Monocyte chemotactic protein MCP-2 activates human basophil and eosinophil leukocytes similar to MCP-3. *J. Immunol.* 154: 4166-4172.
4. Combadiere, C., et al. 1995. Monocyte chemoattractant protein-3 is a functional ligand for CC chemokine receptors 1 and 2B. *J. Biol. Chem.* 270: 29671-29675.
5. Proost, P., et al. 1996. Human monocyte chemotactic proteins-2 and -3: structural and functional comparison with MCP-1. *J. Leukoc. Biol.* 59: 67-74.
6. Dubois, P.M., et al. 1996. Early signal transduction by the receptor to the chemokine monocyte chemotactic protein-1 in a murine T cell hybrid. *J. Immunol.* 156: 1356-1361.

CHROMOSOMAL LOCATION

Genetic locus: CCL13 (human) mapping to 17q11.2.

SOURCE

MCP-4 (E-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping at the C-terminus of MCP-4 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-9656 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

MCP-4 (E-20) is recommended for detection of MCP-4 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for MCP-4 siRNA (h): sc-72122, MCP-4 shRNA Plasmid (h): sc-72122-SH and MCP-4 shRNA (h) Lentiviral Particles: sc-72122-V.


Molecular Weight of MCP-4: 7 kDa.

Positive Controls: MCP-4 (h3): 293T Lysate: sc-110922.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Western blot analysis of human recombinant MCP-4 (A,B). Antibodies tested include MCP-4 (C-18): sc-9656 (A) and MCP-4 (E-20): sc-9656 (B).



MCP-4 (E-20): sc-9656. Western blot analysis of MCP-4 expression in non-transfected (A) and human MCP-4 transfected: sc-110922 (B) 293T whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.


 MONOS
Satisfaction
Guaranteed

Try **MCP-4 (D-2): sc-271124**, our highly recommended monoclonal alternative to MCP-4 (E-20).