

PTH (N-18): sc-9676

BACKGROUND

Parathyroid hormone (PTH), which is also designated parathyrin, is an 84 amino acid single chain peptide that functions to regulate calcium metabolism by raising blood levels of calcium through various mechanisms. PTH stimulates bone formation to increase bone mass and strength in rats and humans. Within the PTH molecule, the essential activity is associated with the first 34 amino acids at the amino terminus of the molecule. The gene encoding PTH maps to human chromosome 11p15.2. Parathyroid hormone-related protein (PTH-rP) is an autocrine factor that is structurally related to PTH yet, unlike PTH, which is synthesized only by the parathyroid cells, PTH-rP is synthesized by several cell types. PTH-rP regulates endochondral bone development and epithelial-mesenchymal interactions during the formation of the mammary glands and teeth. Isolated from the culture medium of a human lung cancer cell line, PTH-rP produces PTH-like effects that are characterized as humoral hypercalcemia of malignancy. The gene encoding PTH-rP maps to human chromosome 12p12.1-p11.2. PTH and PTH-rP are both regulated by vitamin D and steroid hormones and preferentially bind to specific PTH/PTH-rP receptors, then activating adenylate cyclase or PLC β via PKC activation.

CHROMOSOMAL LOCATION

Genetic locus: PTH (human) mapping to 11p15.2; Pth (mouse) mapping to 7 F1.

SOURCE

PTH (N-18) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the N-terminus of PTH of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-9676 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

PTH (N-18) is recommended for detection of precursor and mature PTH of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PTH (N-18) is also recommended for detection of precursor and mature PTH in additional species, including equine, canine and porcine.


Suitable for use as control antibody for PTH siRNA (h): sc-39693, PTH siRNA (m): sc-39694, PTH shRNA Plasmid (h): sc-39693-SH, PTH shRNA Plasmid (m): sc-39694-SH, PTH shRNA (h) Lentiviral Particles: sc-39693-V and PTH shRNA (m) Lentiviral Particles: sc-39694-V.

Molecular Weight of PTH: 9 kDa.


RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941. 3) Immunohistochemistry: use ImmunoCruz™: sc-2053 or ABC: sc-2023 goat IgG Staining Systems.

DATA



PTH (N-18): sc-9676. Mouse thyroid and parathyroid tissue; PTH Green, Keratin 14 Red. Formalin fixed paraffin section with tris retrieval. Kindly provided by A.G. Farr, University of Washington, and M.C. Zúñiga University of California Santa Cruz.



PTH (N-18): sc-9676. Immunoperoxidase staining of formalin fixed, paraffin-embedded human parathyroid gland tissue showing membrane and cytoplasmic staining of glandular cells and extracellular staining of connective tissue.

SELECT PRODUCT CITATIONS

- Sharon, J.L., et al. 2008. The use of N-terminal immobilization of PTH(1-34) on PLGA to enhance bioactivity. *Biomaterials* 29: 3137-3142.
- Dooley, J., et al. 2009. Lessons from thymic epithelial heterogeneity: FoxN1 and tissue-restricted gene expression by extrathymic, endodermally derived epithelium. *J. Immunol.* 183: 5042-5049.
- Kamitani-Kawamoto, A., et al. 2011. MafB interacts with Gcm2 and regulates parathyroid hormone expression and parathyroid development. *J. Bone Miner. Res.* 26: 2463-2472.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

MONOS
Satisfaction
Guaranteed

Try **PTH (H-7): sc-398856** or **PTH (BGN/1F8): sc-80924**, our highly recommended monoclonal alternatives to PTH (N-18).