

G2A (Y-20): sc-9694

BACKGROUND

G2A (for G2 accumulation) is a seven transmembrane G protein-coupled receptor that is up-regulated in response to DNA damage and stress. G2A is predominantly expressed in hematopoietic tissues and in hematopoietic stem cells, and it is more highly detected in pro-B cells, while lower expression is observed in immature B cells and pre-B cells. G2A is expressed throughout T cell maturation, and it is further increased in response to T cell activation. Ectopic expression of a G2A fusion protein in NIH/3T3 fibroblasts induces a cell cycle arrest that is consistent with a block at the G₂/M transition. G2A is also able to attenuate the proliferative effects of Bcr-Abl, a chimeric tyrosine kinase oncogene, suggesting that G2A possesses anti-oncogenic properties. The amino acid sequence of G2A contains a destruction box motif that is consistently observed in cyclins, where it is required for ubiquitination and proteolytic degradation.

REFERENCES

1. Bedi, A., Barber, J.P., Bedi, G.C., el-Deiry, W.S., Sidransky, D., Vala, M.S., Akhtar, A.J., Hilton, J. and Jones, R.J. 1995. BCR-ABL-mediated inhibition of apoptosis with delay of G₂/M transition after DNA damage: a mechanism of resistance to multiple anticancer agents. *Blood* 86:1148-1158.
2. Hochstrasser, M. 1996. Ubiquitin-dependent protein degradation. *Annu. Rev. Genet.* 30: 405-439.
3. Allday, M.J., Inman, G.J., Crawford, D.H. and Farrell, P.J. 1995. DNA damage in human B cells can induce apoptosis, proceeding from G1/S when p53 is transactivation competent and G₂/M when it is transactivation defective. *EMBO J.* 14: 4994-5005.
4. Weng, Z., Fluckiger, A.C., Nisitani, S., Wahl, M.I., Le, L.Q., Hunter, C.A., Fernal, A.A., Le Beau, M.M. and Witte, O.N. 1998. A DNA damage and stress inducible G protein-coupled receptor blocks cells in G₂/M. *Proc. Natl. Acad. Sci. USA* 95: 12334-12339.
5. Shimizu, A., Nishida, J., Ueoka, Y., Kato, K., Hachiya, T., Kuriaki, Y. and Wake, N. 1998. CyclinG contributes to G₂/M arrest of cells in response to DNA damage. *Biochem. Biophys. Res. Commun.* 242: 529-533.
6. Aguda, B.D. 1999. A quantitative analysis of the kinetics of the G₂ DNA damage checkpoint system. *Proc. Natl. Acad. Sci. USA* 96: 11352-11357.

SOURCE

G2A (Y-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping near the C-terminus of G2A of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-9694 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

G2A (Y-20) is recommended for detection of G2A of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for G2A siRNA (m): sc-44371, G2A shRNA Plasmid (m): sc-44371-SH and G2A shRNA (m) Lentiviral Particles: sc-44371-V.

Molecular Weight of G2A: 42-46 kDa.

Positive Controls: WEHI-231 whole cell lysate: sc-2213 or rat skeletal muscle tissue extract: sc-364810.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **G2A (G-5): sc-137112**, our highly recommended monoclonal alternative to G2A (Y-20).