

Neuralized-2 (N-13): sc-98033

BACKGROUND

Neuralized-2, also known as NEURL2, OZZ or Ozz-E3, is a 285 amino acid cytoplasmic protein involved in myofiber differentiation and maturation that is induced during differentiation. Expressed predominantly in skeletal and cardiac muscle, Neuralized-2 localizes to the tips of differentiating myofibers and regulates ubiquitination and degradation of β -catenin at the sarcolemma. Neuralized-2 contains one SOCS box domain, through which it interacts with Elongin B and Elongin C, and one NHR (neuralized homology repeat) domain that is implicated in ubiquitination substrate binding. The gene encoding Neuralized-2 maps to human chromosome 20, which is comprised of approximately 2% of the human genome, contains nearly 63 million bases and encodes over 600 genes. Creutzfeldt-Jakob disease, amyotrophic lateral sclerosis and Alagille syndrome are associated with chromosome 20.

REFERENCES

1. Mead, S., Beck, J., Dickinson, A., Fisher, E.M. and Collinge, J. 2000. Examination of the human prion protein-like gene doppel for genetic susceptibility to sporadic and variant Creutzfeldt-Jakob disease. *Neurosci. Lett.* 290: 117-120.
2. Nastasi, T., Bongiovanni, A., Campos, Y., Mann, L., Toy, J.N., Bostrom, J., Rottier, R., Hahn, C., Conaway, J.W., Harris, A.J. and D'Azzo, A. 2004. Ozz-E3, a muscle-specific ubiquitin ligase, regulates β -catenin degradation during myogenesis. *Dev. Cell* 6: 269-282.
3. Nishimura, A.L., Mitne-Neto, M., Silva, H.C., Oliveira, J.R., Vainzof, M. and Zatz, M. 2004. A novel locus for late onset amyotrophic lateral sclerosis/motor neurone disease variant at 20q13. *J. Med. Genet.* 41: 315-320.
4. Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 608597. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Song, R., Koo, B.K., Yoon, K.J., Yoon, M.J., Yoo, K.W., Kim, H.T., Oh, H.J., Kim, Y.Y., Han, J.K., Kim, C.H. and Kong, Y.Y. 2006. Neuralized-2 regulates a Notch ligand in cooperation with Mind bomb-1. *J. Biol. Chem.* 281: 36391-36400.

CHROMOSOMAL LOCATION

Genetic locus: NEURL2 (human) mapping to 20q13.12.

SOURCE

Neuralized-2 (N-13) is an affinity purified goat polyclonal antibody raised against a peptide mapping within a cytoplasmic domain of Neuralized-2 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-98033 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

Neuralized-2 (N-13) is recommended for detection of Neuralized-2 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Neuralized-2 siRNA (h): sc-75905, Neuralized-2 shRNA Plasmid (h): sc-75905-SH and Neuralized-2 shRNA (h) Lentiviral Particles: sc-75905-V.


Molecular Weight of Neuralized-2: 31 kDa.

Positive Controls: LADMAC whole cell lysate: sc-364189.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



Neuralized-2 (N-13): sc-98033. Western blot analysis of Neuralized-2 expression in KNRK (A) and LADMAC (B) whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **Neuralized-2 (A-11): sc-398195**, our highly recommended monoclonal alternative to Neuralized-2 (N-13).