

HIRIP3 (M-43): sc-98402

BACKGROUND

The HIRIP3 (HIRA interacting protein 3) locus encodes for a 556 amino acid protein that directly interacts with the HIRA histone chaperone. It also interacts weakly with core histones, Histone H2B and Histone H3. HIRIP3 is a heavily phosphorylated nuclear protein and it is found throughout the cell cycle. It is phosphorylated by casein kinase II. HIRIP3 may play a role in chromatin function and histone metabolism. A region (approximately 60 amino acids in length) at the C-terminus of HIRIP3 is highly conserved among vertebrates and it contains residues that are invariantly charged, polar and hydrophobic. Two isoforms of HIRIP3 exist due to alternative splicing. Isoform 1 is predominantly expressed in skeletal muscles and isoform 2 is expressed in the liver and the heart. Human HIRA homologs are thought to be responsible for the DiGeorge syndrome and related developmental disorders.

REFERENCES

- Lorain, S., Quivy, J.P., Monier-Gavelle, F., Scamps, C., Lécluse, Y., Almouzni, G. and Lipinski, M. 1998. Core histones and HIRIP3, a novel histone-binding protein, directly interact with WD repeat protein HIRA. *Mol. Cell. Biol.* 18: 5546-5556.
- Magnaghi, P., Roberts, C., Lorain, S., Lipinski, M. and Scambler, P.J. 1998. HIRA, a mammalian homologue of *Saccharomyces cerevisiae* transcriptional co-repressors, interacts with Pax3. *Nat. Genet.* 20: 74-77.
- Online Mendelian Inheritance in Man, OMIM™. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 603365. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Lorain, S., Lécluse, Y., Scamps, C., Mattéi, M.G. and Lipinski, M. 2001. Identification of human and mouse HIRA-interacting protein-5 (HIRIP5), two mammalian representatives in a family of phylogenetically conserved proteins with a role in the biogenesis of Fe/S proteins. *Biochim. Biophys. Acta* 1517: 376-383.
- Ahmad, A., Kikuchi, H., Takami, Y. and Nakayama, T. 2005. Different roles of N-terminal and C-terminal halves of HIRA in transcription regulation of cell cycle-related genes that contribute to control of vertebrate cell growth. *J. Biol. Chem.* 280: 32090-32100.
- Frith, M.C., Ponjavic, J., Fredman, D., Kai, C., Kawai, J., Carninci, P., Hayashizaki, Y., Hayshizaki, Y. and Sandelin, A. 2006. Evolutionary turnover of mammalian transcription start sites. *Genome Res.* 16: 713-722.
- Luk, E., Vu, N.D., Patteson, K., Mizuguchi, G., Wu, W.H., Ranjan, A., Backus, J., Sen, S., Lewis, M., Bai, Y. and Wu, C. 2007. Chz1, a nuclear chaperone for Histone H2AZ. *Mol. Cell* 25: 357-368.
- Assrir, N., Filhol, O., Galisson, F. and Lipinski, M. 2007. HIRIP3 is a nuclear phosphoprotein interacting with and phosphorylated by the serine-threonine kinase CK2. *Biol. Chem.* 388: 391-398.
- SWISS-PROT/TrEMBL (Q9BW71). World Wide Web URL: <http://www.expasy.ch/sprot/sprot-top.html>

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: Hrip3 (mouse) mapping to 7 F3.

SOURCE

HIRIP3 (M-43) is a rabbit polyclonal antibody raised against amino acids 46-88 mapping near the N-terminus of HIRIP3 of mouse origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

HIRIP3 (M-43) is recommended for detection of HIRIP3 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for HIRIP3 siRNA (m): sc-145972, HIRIP3 shRNA Plasmid (m): sc-145972-SH and HIRIP3 shRNA (m) Lentiviral Particles: sc-145972-V.

Molecular Weight (predicted) of HIRIP3: 62 kDa.

Molecular Weight (observed) of HIRIP3: 90 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.


 MONOS
Satisfaction
Guaranteed

Try **HIRIP3 (E-9): sc-514646**, our highly recommended monoclonal alternative to HIRIP3 (M-43).