

PCP-4 (FL-62): sc-98549

BACKGROUND

PCP-4 (Purkinje cell protein 4), also known as PEP-19, is a calmodulin (CaM) regulatory protein that is highly expressed in neuronal cells. Through its IQ motif, PCP-4 mediates both the calcium-dependent binding properties of CaM and the rates of association and dissociation of calcium from the C-terminal domain of CaM. The IQ motif contains a serine residue which can be phosphorylated by all isoforms of protein kinase C (PKC). PCP-4 is implicated in uterine leiomyomas as well as in neurodegenerative disorders such as Alzheimer's disease and Huntington's disease. Additionally, overexpression of PCP-4 is thought to play a role in cerebellar hypoplasia, a key feature of Down syndrome.

REFERENCES

- Erhardt, J.A., Legos, J.J., Johanson, R.A., Slemmon, J.R. and Wang, X. 2000. Expression of PEP-19 inhibits apoptosis in PC12 cells. *Neuroreport* 11: 3719-3723.
- Slemmon, J.R., Feng, B. and Erhardt, J.A. 2000. Small proteins that modulate calmodulin-dependent signal transduction: effects of PEP-19, neuro-modulin, and neurogranin on enzyme activation and cellular homeostasis. *Mol. Neurobiol.* 22: 99-113.
- Putkey, J.A., Kleerekoper, Q., Gaertner, T.R. and Waxham, M.N. 2003. A new role for IQ motif proteins in regulating calmodulin function. *J. Biol. Chem.* 278: 49667-49670.
- Kanamori, T., Takakura, K., Mandai, M., Kariya, M., Fukuhara, K., Kusakari, T., Momma, C., Shime, H., Yagi, H., Konishi, M., Suzuki, A., Matsumura, N., Nanbu, K., Fujita, J. and Fujii, S. 2003. PEP-19 overexpression in human uterine leiomyoma. *Mol. Hum. Reprod.* 9: 709-717.
- Skold, K., Svensson, M., Nilsson, A., Zhang, X., Nydahl, K., Caprioli, R.M., Svenningsson, P. and Andr n, P.E. 2006. Decreased striatal levels of PEP-19 following MPTP lesion in the mouse. *J. Proteome Res.* 5: 262-269.
- Simons, M.J. and Pellionisz, A.J. 2006. Genomics, morphogenesis and biophysics: triangulation of Purkinje cell development. *Cerebellum* 5: 27-35.
- Dickerson, J.B., Morgan, M.A., Mishra, A., Slaughter, C.A., Morgan, J.I. and Zheng, J. 2006. The influence of phosphorylation on the activity and structure of the neuronal IQ motif protein, PEP-19. *Brain Res.* 1092: 16-27.
- Sala, A., Scaturro, M., Proia, P., Schiera, G., Balistreri, E., Aflalo-Rattenbach, R., Creau, N. and Di Liegro, I. 2007. Cloning of a rat-specific long PCP-4/PEP-19 isoform. *Int. J. Mol. Med.* 19: 501-509.
- Xiao, J., Wu, Y., Chen, R., Lin, Y., Wu, L., Tian, W. and Liu, L. 2008. Expression of PCP-4 gene during osteogenic differentiation of bone marrow mesenchymal stem cells *in vitro*. *Mol. Cell. Biochem.* 309: 143-150.

CHROMOSOMAL LOCATION

Genetic locus: PCP4 (human) mapping to 21q22.2; Pcp4 (mouse) mapping to 16 C4.

SOURCE

PCP-4 (FL-62) is a rabbit polyclonal antibody raised against amino acids 1-62 representing full length PCP-4 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

PCP-4 (FL-62) is recommended for detection of PCP-4 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

PCP-4 (FL-62) is also recommended for detection of PCP-4 in additional species, including canine, bovine and avian.

Suitable for use as control antibody for PCP-4 siRNA (h): sc-76091, PCP-4 siRNA (m): sc-76092, PCP-4 shRNA Plasmid (h): sc-76091-SH, PCP-4 shRNA Plasmid (m): sc-76092-SH, PCP-4 shRNA (h) Lentiviral Particles: sc-76091-V and PCP-4 shRNA (m) Lentiviral Particles: sc-76092-V.

Molecular Weight of PCP-4: 7 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



MONOS
Satisfaction
Guaranteed

Try **PCP-4 (1E3): sc-293258**, our highly recommended monoclonal alternative to PCP-4 (FL-62).